

BOSCAROL

Emergency Systems

MILITARY MEDICAL EQUIPMENTS



MILITARY LINE



MEDCAL BAG

- Compact but really capacious emergency bag
- 4 independent pockets, 3 of which with zip
- Made in tearproof impermeable camouflage material
- 2 rings for transport on shoulders

ART. COD

MIL 10000

MEDICAL BAG



MEDCAL BAG GREEN

- Compact but really capacious emergency bag
- 4 independent pockets, 3 Made in tearproof impermeable green material
- 2 rings for transport on shoulders
- 1 loop for belt
- 1 reflective strips

ART. COD

MIL 10004

MEDICAL BAG GREEN




RESCUE BAG

- Compact and practical rescue bag
- Easy to carry either by hand, on shoulder or as a backpack
- Made of shockproof impermeable camouflage material
- 2 windows for identity cards/logos insertion
- Lower part with inner basin suitable for the oxygen cylinder
- Upper part with elastics
- Dimensions: about 47x26x27 cm

ART. COD

MIL 10010

MILITARY RESCUE BAG



TRAUMABAG 10 MILITARY

- Functional, compact and capacious emergency bag.
- Transportable by hand, on shoulder or as a backpack
- 2 external pockets, frontal and rear equipped with elastic
- 2 lateral external pockets
- 4 inner modular spaces
- Removable coloured pouch with transparent frontal
- Made of tearproof impermeable camouflage
- Entirely removable for easy washing
- Dimensions: 45x30x31h cm
- Weight: 2,1 Kg

ART. COD

MIL 10020

MILITARY TRAUMA BAG 10



TRAUMABAG 20

- Functional, compact and capacious emergency bag.
- Transportable by hand, on shoulder or as a backpack
- 2 external pockets, anterior one equipped with elastics
- 2 lateral external pockets
- 5 internal modular spaces
- Removable coloured pouch with transparent frontal
- Made of tearproof impermeable camouflage
- Entirely removable for easy washing
- Dimensions: about 53x33x34h cm
- Weight: 2,3 Kg

ART. COD

MIL 10026

MILITARY TRAUMA BAG 20



TRAUMABAG 30

- Functional, compact and capacious emergency bag.
- Transportable by hand, on shoulder or as a backpack
- 3 external pockets, 2 of which frontal and equipped with elastics
- 2 lateral external pockets
- 5 internal modular spaces
- Removable coloured pouch with transparent frontal
- Made of tearproof impermeable camouflage
- Entirely removable for easy washing
- Dimensions: about 63x37x31h cm
- Weight: 2,7 Kg

ART. COD

MIL 10032

TRAUMA BAG 30

MILITARY LINE



CE



CE



CE

MILITARY PC RUCKSACK

- Compact emergency rucksack
- Divided into 2 parts: the lower part is rigid and protected, the upper part has a rigid perimeter structure
- Complete with pockets, compartments and elastics for medical material
- Internal book-system pockets and elastics for further material
- Made in tear-proof impermeable camouflage material
- 1 frontal window allowing identity cards/logos insertion
- Dimensions: about 45x24x48h cm
- Weight: 2,4Kg

EMERGENCY RUCKSACK VIPER

- Rucksack for advanced emergency
- Compact and practical
- Complete with 9 coloured pouches with transparent frontal and Velcro strips
- Internal Velcro strips for further pouches fixing.
- 2 frontal pockets with zip
- 3 external perimeter pockets
- Elastics all along bag perimeter suitable for further medical material fixing
- Central dividing bulkhead with zip, complete with elastics
- Made in tear-proof impermeable camouflage
- Dimensions: about 42x23x51 cm
- Weight: 2,8 Kg

ATLS RUCKSACK

- A TLS emergency rucksack
- 1 wide frontal pocket
- Rigid red shockproof inner basin with central oxygen cylinder holder.
- 1 zipped compartment with Velcro strips
- Wide internal pocket
- Made in tear-proof impermeable camouflage
- P added shoulder-straps
- Dimensions: about 32x22x52h cm

ART. COD

MIL 10220

MILITARY PC RUCKSACK

ART. COD

MIL 10200

VIPER EMERGENCY RUCKSACK

ART. COD

MIL 10240

ATLS RESCUE RUCKSACK



CE



CE



CE

MODULE

- Great capacity emergency rucksack
- 2 frontal external pocket
- 5 internal removable coloured pouches
- 4 coloured pouches with transparent frontal
- Internal Velcro strips for the fixing of further pouches
- Realized in tearless impermeable camouflage
- Dimensions: 39x36x59h cm
- Weight: 3,6 Kg

RESCUE RUCKSACK 02

- Made of tear-proof impermeable camouflage
- Big frontal pocket
- 1 zip fastened window allowing insertion of identity cards/logos
- Inner red, rigid, shockproof basin
- Lateral oxygen cylinder holder
- Shoulder-straps for transport as a rucksack
- Dimensions: about 27x20x47h cm

RESCUE RUCKSACK

- Made of tear-proof impermeable camouflage
- Big frontal pocket
- 1 zip fastened window allowing insertion of identity cards/logos
- Inner red, rigid, shockproof basin
- 4 coloured pouches with transparent frontal and handle
- Shoulder-straps for transport as a rucksack
- Dimensions: about 27x20x47h cm
- Weight: 1,8 Kg

ART. COD

MIL 10230

EMERGENCY RUCKSACK
MODULE CAMOUFLAGE

ART. COD

MIL 10212

RESCUE RUCKSACK 02

ART. COD

MIL 10210

RESCUE RUCKSACK

MILITARY LINE



MASIER RUCKSACK

- 2 frontal pockets equipped with transparent, zips and elastics
- Inner red, rigid, shockproof basin with central support for oxygen cylinder
- 1 compartment with zip and Velcro-strips
- 1 inner wide transparent pocket
- 4 removable coloured pouches with transparent frontal
- Made in tear-proof impermeable camouflage material
- Padded shoulder-straps
- Dimensions: about 32x22x52h cm

ART. COD

MIL 10260

MILITARY MASTER EMERGENCY RUCKSACK



RUCKSACK MAXI

- Rucksack for advanced rescue
- Compact and practical
- 2 frontal pockets equipped as drugs bag holder (the first one) and for infusions (the second one)
- 2 external pockets: 1 arranged with 4 coloured pouches with transparent frontal for diagnostics, the other one for the reanimation bag
- Inner red, rigid, shockproof basin complete with 3 coloured pouches with transparent frontal
- 4 coloured transparent pouches with transparent frontal and Velcro fixed in inside
- Inner Velcro strips for the fixing of further pouches
- Made in tear-proof impermeable

ART. COD

MIL 10250

MAXI RESCUE RUCKSACK



MEDBAG

- Suitable for use by Medics & Paramedics in emergency & first aid scenarios.
- Lightweight, yet with aluminium case practicality.
- Realized with a semirigid basin covered with tearproof impermeable camouflage
- Washable tearproof thermo-shaped internal with modular compartments
- Arranged for both 1 and 2 litres oxygen cylinder storage.
- Upper part complete with 4 removable coloured pouches with transparent frontal
- 2 dividing wall with Velcro for the fixing of further optional pouches
- Dimensions: 53x38x23 cm
- Weight: 3,8 Kg

ART. COD

MIL 10400

MEDIBAG



INTUBATION BAG

- Bag for the intubation
- Interior logic structure arranged with elastics for the intubation set fastening
- Made in tear-proof impermeable camouflage material
- 2 Velcro-strips for the fixing on rucksacks or rescue bags
- Dimensions: about 32x7,5x22h cm

ART. COD

MIL 10520

INTUBATION BAG



MAXIDRUGS BAG WITH ELASTICS

- Book-opening drugs bag
- Transport handle
- 2 removable plates with elastics for a total of 106 different sized phials
- Interior padded division with transparent list holder for phials expiry date
- Made in tear-proof impermeable camouflage material
- 2 Velcro-strips for the fixing on bags or rucksacks.
- Dimensions: about 46x9x34 cm

ART. COD

MIL 10540

MAXI DRUGS BAG WITH ELASTICS



MAXIDRUGS CASE WITH RAIL

- Drugs bag with book-opening
- 2 removable plates able to contain a total of 104 phials
- 8 RAIL ampoules holders for 13 positions
- Interior padded divider with pockets for needles and further medical material
- Made in tear-proof impermeable camouflage material
- 2 Velcro-strips for the fixing on bags or emergency rucksacks

ART. COD

MIL 10544

MAXI DRUGS BAG WITH RAIL

MILITARY LINE



CE

DRUGSBAG WITH ELASTICS

- Drug bag with book-opening
- 2 removable plates with elastics able to contain a total of 72 phials
- Interior padded divider with pockets for needles and further medical material
- Made in tear-proof impermeable camouflage material
- 2 Velcro-strips for the fixing on bags or emergency rucksacks

ART. COD

MIL 10530

DRUGS BAG WITH ELASTICS



CE

RAIL DRUGS BAG

- Drug bag with book-opening
- 2 removable plates with rail able to contain up to 60 phials
- Interior padded divider with pockets for needles and further medical material
- Made in tear-proof impermeable camouflage material
- 2 Velcro-strips for the fixing on bags or emergency rucksacks
- Dimensions: about 32x7,5x22,5 cm

ART. COD

MIL 10534

RAIL DRUGS BAG



CE

PROF® VACUUM MATTRESS

- New orthopaedic PROF® vacuum mattress with 3 independent chambers
- Realised in compliance with EN 1789/1999 normative
- CE certified in accordance with 1865/1999 norm
- Perfect and absolute immobilization of injured person
- High-resistance cloth
- Protective cloth made in tear-proof impermeable camouflage material
- Compact suction valve with 1/4 on/off
- 3 handles each side obtained from the protective cloth itself
- 2 head-handles and 2 foot-handles.
- 3 safety belts
- 100% X-ray transparent.
- CAT scan compatible
- Dimensions: 210x100 cm

ART. COD

MIL 10600

PROF® VACUUM MATTRESS



CE

2 ANATOMIC VACUUM SPLINTS

- Set of 2 anatomic vacuum splints suitable for leg, arm and ankle
- Absolute, quick and painless immobilization of fractured limbs
- High resistant cloth
- External, compact suction valve with 1/4 on/off
- Protection pockets made of tearproof, impermeable camouflage cloth with fixing system and rapid hook/unhook moulding
- They ensure a perfect immobilization without compressing, ensuring physiological blood circulation
- 100% X-ray translucent
- CAT scan compatible
- Complete with transport bag
- Arm/ankle dimensions: about 68x38x53

ART. COD

MIL 10610

2 ANATOMIC VACUUM SPLINTS SET WITH BAG



CE

TRANSFER SHEET

- Transfer sheet made of extremely resistant, impermeable, washable and disinfectable cloth.
- 8 handles allowing an easy transport of patient.
- Camouflage cloth
- Dynamic resistance at 680 Kg/cm²
- Dimensions: 2000x700 mm

ART. COD

MIL 10620

TRANSFER SHEET



CE

AS PANTS

- Anti-shock pants for patients experiencing hypovolemic shock.
- Easy application.
- 3 independent chambers.
- 3 Ventil valves with 1/4 on/off.
- 3 security valves calibrated to 105 mm Hg.
- 3 monitoring manometers with protection.
- Even pressure.
- Large groin-opening.
- Internal tubes.
- 1 single connection for the pump and the 3 chambers.
- Chambers protection made in tearproof camouflage cloth.
- 1 inflation pump.
- 1 transport case in military camouflage cloth.

ART. COD

MIL 10630

AS PANTS

MILITARY LINE



CE

TRACTIONSPLINT SET

- Traction and immobilization system for thigh- bone fractures
- It alleviates pain and avoids further muscular or nervous damages around the fracture
- Aluminium structure adjustable in length, easy to use
- Belts with Velcro closing
- Soft ankle for ankle fixing
- Automatic adjustable traction mechanism, allowing broken thighbone alignment in its natural position
- Complete set featuring adult and pediatric traction splint with transport bag
- Dimensions: about 107x38 cm
- Weight: 4,7Kg

ART. COD

IMM 123060

TRACTION SPLINT SET



CE

PNEUMATIC SPLINT SET -

- Set of 6 pneumatic splints.
- Short and long leg, short and long arm, ankle, hand.
- Easy utilization through a zipper.
- Simple inflation valve.
- Realized in PVC material, resistant from - 25° to + 45°
- Supplied with pump.
- Practical transport bag.
- Dimensions set: about 38x5x39h. cm
- Weight: 1,7 Kg.

ART. COD

IMM 120270

PNEUMATIC SPLINT SET - 6 PIECES



CE

OB2012 MILITARY

OB2024 MILITARY

- Medical Suction Unit OB 2012 powerful, easy and compact.
- Silent and powerful double piston motor of more than 30 lt./min.
- Adjustable linear suction, with control-manometer from 0 to > 800 mbar.
- Autonomy of about 45 Minutes.
- New box made in ABS/UL.
- Autoclavable 1000 ml-jar for the secretion collection with overflow valve and filter directly integrated in the cover.
- Tear and water-proof transport bag made in black Tyropol
- Cable for the connection to the 12Volts of the vehicle.
- Dimensions: 35x12x24 cm.
- Weight: about 4,5 Kg

ART. COD

MIL 10800

OB 2012 MILITARY MEDICAL SUCTION UNIT

ART. COD

MIL 10804

OB 2024 MILITARY MEDICAL SUCTION UNIT



CE

OB 1000 MEDICAL SUCTION

- Medical suction unit OB 1000, compact and powerful.
- It operates directly through the 12 Volt C.C of the vehicle and through the internal accumulator with an autonomy of over 45 minutes.
- Motor of over 20 lt./min capacity
- Suction superior to 750mbar
- Micrometrical vacuum regulation through vacuum control gauge
- Box realised in ABS/UL
- 1000 ml autoclavable jar with filter directly integrated in the closing
- Supplied complete of tearproof impermeable camouflage transport bag
- Connecting cable to 12Vcc of the vehicle
- Dimensions: cm 31x28x8h
- Weight: 3,5 Kg

ART. COD

MIL 10820

OB 1000 MILITARY MEDICAL SUCTION UNIT



CE

MILITARY OB MINIVAC

- OB MINIVAC, compact medical military suction unit:
- It functions directly through the vehicle 12V C.C.
- Over 45 minutes autonomy
- Adjustable suction: over 700 mbar
- Control vacuum gauge
- Box realised in ABS/UL
- 500 ml autoclavable jar for the secretion collection with external filter
- 12 V connecting cable
- Serial battery charger
- Green transport bag
- Dimensions: 28x7x20 cm
- Weight: 2,2 kg

ART. COD

MIL 10840

MILITARY OB MINIVAC



CE

OXBAG

- OXIBAG padded cylinder-holder bag with handle and fixing system to the stretcher.
- Thanks to its particular shape and padding it can safely contain, secure and transport oxygen cylinder and pressure reducer.
- It can be safely fastened to the stretcher according to standard EN 1789.
- Made of high resistant tearproof impermeable camouflage cloth and may contain a 2 litres oxygen cylinder

ART. COD

MIL 10500

CAMOUFLAGE OXIBAG

MILITARY LINE



CE

BPACK MILITARY

- B PPACK two rucksacks in one
- Capacious multipurpose rescue rucksack BPACK.
- Complete with 6 coloured removable pouches fixed with Velcro
- External superior pouch
- Independent mini-rucksack fixed on the front and completely removable.
- 3 perimetric external pouches
- Anatomic back
- Adjustable shoulder-straps
- Realized in tearproof camouflage material



CE

MILITARY MULTIPACK

- New emergency rucksack, realized in collaboration with the Helicopter Rescue service.
- Rigid structure rucksack with protected perimetric borders and 2 transport handles
- 7 removable rigid pouches with transparent frontal and zipper.
- 1 long lateral pouch suitable for catheters storage
- Wide frontal pocket with zipper
- Stuffed shoulder straps contained in its proper posterior pocket.
- Realized in tearproof impermeable camouflage cloth
- Dimensions: cm 35x25x68h
- Weight: about 5,1 kg



CE

WARMING INFUSION BAG

- Special infusion warming allowing heating and maintenance of 2 infusions temperature up to 37°
- Inner cloth with carbon fibres working at 12 Volt.
- In case of emergency activities, an inner pocket for the possible insertion of an instantaneous or reusable warm tablet is present.
- Realized in blue Cordura material. Inner covering made in insulating thermo-refractive material
- Lateral transparent window suitable for examination.
- A posterior pocket suitable for flow control devices storage.
- 12 Volt cable with lighter plug
- Dimensions: cm 23x13x28h
- Complete of THERMIC PROTECTION for IV ADMINISTRATION cm 150

ART. COD

MIL 10280

BPACK MILITARY

ART. COD

MIL 10290

MILITARY MULTIPACK
EMERGENCY RUCKSACK

ART. COD

FCT 2314

WARMING INFUSION BAG WITH
THERMIC PROTECTION



CE

BURN KIT 'N POUCH

- Burn kit, burn treatment set in equipped tearproof impermeable material with 2 coloured pouches with handle, complete with:
- N o. 01 VESTATEX cover 70x80 cm, 1300 g. Hydrogel
- N o. 01 VESTASAN cloth 35x50, 360 g. Hydrogel
- N o. 02 VESTASAN gauzes 10x10 cm with 40 g Hydrogel
- N o. 01 VESTASAN gauze 10x15 cm with 100 gm Hydrogel
- N o. 02 VESTAGEL 25 ml Hydrogel sachets • N o. 01 100 ml Hydrogel Vesta bottle
- N o. 02 pairs of sterile gloves
- N o. 02 elastic medication bandages
- N o. 01 emergency scissor, plastic handle
- N o. 02 instant ice packs
- N o. 02 not injectable sodium chloride 250 ml bottles

ART. COD

VAL 40510

BURN KIT

ART. COD

VAL 44052

B.I.G. SET ADULT



CE

B.I.G. SET ADULT

- B .I.G. device allowing easy INTRABONE access, fast and safe when venous access is not possible.
- Set in pouch complete with sterile gloves, lodopovidone, plaster, gauzes.



CE

ROLLSPLINT

- Adjustable light splint suitable for fractures, washable and disinfectable.
- Covered with atoxic material.
- X ray translucent
- Dimensions: cm 92x11

ART. COD

IMM 123020

ROLLSPLINT

MILITARY LINE



CE

ISOTHERMIC BLANKET

- Protection and rescue isothermic blanket
- Prevents patients' hypothermia whether traumatized, burned or post-operated
- Double function: warming and refreshing
- Silver side, if positioned toward the body, isolates and reflects, preserving about 80% of body warmth. Thus protects patient from hypothermia, coldness, freezing and from shock.
- Golden side: placed toward patient body and silver one, reflects sun rays maintaining body freshness.

ART. COD

FCT 3000

ISOTHERMIC BLANKET



CE

WARMING SHEET

- Secure low tension warming sheet 12 Volt AC for therapeutic use.
- Polyester external cloth, internal cloth in spun carbonic fibre achieving the desired temperature in 5 minutes.
- Supplied with 12 V cable
- Warming sheet dimensions:

ART. COD

FCT 2220

WARMING SHEET



CE

EMERGENCY MEDICAL BLANKET

- New emergency blanket
- AUTOCALVABLE
- ANTIBACTERIAL
- TRANSPIRING
- WATER REPELLANT
- WINDPROOF
- THERMIC ISOLATION
- MICROCLIMATE MAINTENANCE
- SOFT, LIGHT, RESISTANT
- Realized in continue filament (tri-cloth)
- Its characteristics made this a unique blanket and renders emergency activity safe and economic.

ART. COD

FCT 2420

EMERGENCY MEDICAL BLANKET



CE

PRESSURE INFUSOR

- Pressure infusor for medicaments contained in plastic material bags. It is complete with manometer for infusion scanning and monoball with deflating valve for ending dosage. Suitable for 500 and 1000 ml. infusion.

ART. COD

DIA 100390

PRESSURE INFUSOR WITH OB MANOMETER



CE

TRIANGULAR SHEET VACUUMPACKED

- Immobilization triangular sheet
- Manufactured in cloth
- Dimensions cm 90x90x129
- 2 safety pins
- Singularly vacuum-packed

ART. COD

MED 96760

TRIANGULAR SHEET VACUUM-PACKED



CE

EXPLORER BAUE

- Rescue case EXPLORER, impermeable and watertight, water resistant up to 10 m. depth
- Neoprene seal and exclusive pressurization valve.
- Made of shockproof resins resisting to corrosive agents
- Absolute protection for the material included.
- Ergonomic handle
- Sturdy click closing arranged for padlock or lead closings.
- 2 trunk types
- Small trunk dimensions: 649x507x354 mm
- Big trunk dimensions: 834x552x330 mm

ART. COD

PC 190140

EXPLORER PC TRUNK S

PC 190144

EXPLORER PC TRUNK B



OUR ENTERPRISES, OUR COMMUNICATION



Internet
www.boscarol.it



General Catalogue on Cd



Four monthly offers
of our products



Means:

1 T.Q.S.
Total Quality Satisfaction
(Quality without

- 2 • Research
- Developme
- nt
- Innovation
- Rationality
- Reliability
- Production

3 Information
and customer
assistance with qualified and prepared

4 International presence

5 Furthermore, we attend as expositor the most important exhibitions, manifestations, presentations and organize dedicated Work Shop in collaboration with our official retailers and agents in order to better show our products range.

Quality means being able to anticipate customer's needs and transfer these needs into a useful and conform product, as well as to create a system able to reduce prices at the lowest possible level so to assure an affair for the customer and yet a profit for the Company. (W. Edwards Deming)

WE OPERATE
IN THE SECTOR OF PREVENTION
TO SUPPORT PEOPLE
WHO EVERY DAY REQUIRE ASSISTANCE AT
WORK:

WE OPERATE
FOR PEOPLE
WHO URGENTLY NEED HELP,
FOR PEOPLE
WHO TRUST OTHER PEOPLE
TO TAKE CARE
OF THEIR SAFETY



HOW TO REACH US AND CONTACT US



39100 BOLZANO - Via G. Di Vittorio, 33 - C.P. 639
 Tel. ++39 0471 932893 - Fax ++39 0471 932331 - Export Dep. ++39 0471 934748 E-mail: info@boscarol.it - www.boscarol.it



BY MAIL
 Use the specific order form, fill in clearly and

BOSCAROL EMERGENCY SYSTEMS S.R.L.
 Via G. Di Vittorio, 33 - P. Box 639 - 39100 BOLZANO



BY FAX
 Fastest and safest way to send us your order. In order to speed up your order please fill in the order form in all its parts

0471 932331




INTERNET
info@boscarol.it



TELEPHONE
 It is possible to receive further information, goods availability and delivery time and worthy advice upon request.

0471 932893

OPENING TIMES
 From Monday to Friday 8.30 am - 12.30 am/2.30 pm - 6.30 pm




FREEWAY EXIT BOLZANO SUD
 FIRST TRAFFIC LIGHT RIGHT
 FIRST ROUNDABOUT RIGHT
 FIRST TRAFFIC LIGHT RIGHT

OUR CERTIFICATIONS



RUSSIAN CONFEDERATION CERTIFICATIONS





NEW ZEALAND POSTAGE
STAMP DEDICATED
TO BOSCAROL SRL



**NATO CODE ASSIGNED
7175P**

We thank the following for the collaboration:

BOSCAROL PER LE SUE SPEDIZIONI SI AFFIDA AL CORRIERE EXECUTIVE



EXECUTIVE
GROUP
INTERNATIONAL

GIÀ CONSEGNATO

Member of  *General E. Parrot*

Sede EXECUTIVE di Bolzano
Via Galvani, 22 - 39100 Bolzano - Italy
Tel. 0471 545000 - Fax 0471 545018
e-mail bolzano@executivegroup.com

Numero Verde
800 33 1393



foto
Gianni
DePolo

340 9245770

LA
Protezione Civile
I T A L I A N A

Direzione generale, Redazione e Amministrazione:
20142 Milano - Viale Faenza, 26/5
Tel. 02 8135018/8136669 - Fax 02 8134925

www.laprotezionecivile.com
e-mail: info@laprotezionecivile.com

Menù di informazione e studi per le componenti del Servizio Nazionale della Protezione Civile



Professionalità
Tecnologia
Qualità


GRAFICA Benvenuto

AREAGRAFICA
PROGETTAZIONE GRAFICA
SVILUPPO IMMAGINE COORDINATA
SOLUZIONI DI STAMPA

AREASTAMPA
STAMPATI COMMERCIALI
RIVISTE, DEPLIANTS, CATALOGHI
POSTERS, CARTELLI VETRINA

Via Varesche, 58
37010 Albarè di Costermano (VR)
Tel.: 045 620 0668 - Fax: 045 720 0433
e-mail: graficabenvenuto@tin.it

THANK YOU PARTICULARLY TO THE PERSON IN CHARGE OF 118, AIR
RESCUE, AND WHITE CROSS OF AUTONOMY
PROVINCE OF BOLZANO (ITALY)



Visit our website
www.boscarol.it
and search for authorized retailers in your



Via G. Di Vittorio, 33
39100 BOLZANO - ITALY
Tel. +39 0471 932893
Fax +39 0471 932331
Export Dep. +39 0471 934748
E-mail: info@boscarol.it
www.boscarol.it

We reserve the right to vary technical specifications without notice.
Accessories are not included. Their purpose is merely descriptive