

V409

Fingertip pulse oximeter

Technology that doesn't miss a beat



The **DARAY** V409 fingertip pulse oximeter. Used to measure oxygen saturation level and pulse rate via the patient's finger

- Simple and comfortable to use
- Clear, wide-angle OLED display
- Rotating display
- Real-time display of SpO2, pulse rate graphical display and PLETH
- Accuracy: SpO2 $\pm 2\%$ and pulse rate ± 2 BPM
- Comes complete with lanyard and batteries
- Extremely lightweight (less than 50g inc. batteries)
- Powered by 2 x AAA alkaline or rechargeable batteries
- Autopower-off after 8 seconds
- ~ 100 hours battery life
- Low-battery indicator

ABOUT DARAY

DARAY have been established in the medical, dental and veterinary markets since 1968 building exceptional and durable products.

We shall continue to innovate; aiming to excel in our endeavours to deliver reliable solutions.



CE
0120

DARAY[®]
Medical

www.daray.com

Tel: 0800 804 8384 Email: info@daray.co.uk