

V501

Fingertip pulse oximeter

Technology that doesn't miss a beat



The **DARAY** V501 fingertip pulse oximeter is used to measure oxygen saturation levels and pulse rate via the patient's finger.

- Simple and comfortable to use
- Graphical pulse rate display
- 50 hours battery life
- Low-battery indicator
- Auto rotation display
- Accuracy $\pm 2\%$ for SpO₂ and ± 2 BP for pulse
- Clear, wide-angle LCD display
- Auto power-off after 8 seconds
- Easy-to-change AAA batteries
- Comes with lanyard, and batteries
- Real-time display of SpO₂
- Extremely lightweight

ABOUT DARAY

DARAY have been established in the medical, dental and veterinary markets since 1968 building exceptional and durable products.

We shall continue to innovate; aiming to excel in our endeavours to deliver reliable solutions.



CE
0120

DARAY[®]
Medical

www.daray.com

Tel: 0800 804 8384 Email: info@daray.co.uk