



ARKY cabinets are the perfect solution for the safe and professional storage of an Automated External Defibrillator

ARKY was founded by Defibrion in 2008. Defibrion is the exclusive distributor of HeartSine AED's in the Benelux countries and France. Since the beginning, ARKY has successfully expanded its business to more than 20 countries worldwide.



POLYSTYRENE OUTDOOR CABINET

High quality

- Made of high impact polystyrene
- Protects the AED against all types of weather
- 100% dust- and waterproof
- Well sealed with brushed stainless steel hinges and closures

Dimensions

- Inside dimensions:
28 x 32 x 18 cm
- Outside dimensions:
34 x 40 x 20 cm

Options

- Acoustic alarm
- Company branding
- Including breakaway seal
- Heating element (24V)



Product no.	Description	Aluminum Inlay	Alarm	Lockable (with breakaway seal/ padlock)	Heating element (24V)
182	Weather resistant plastic cabinet (without heating)	✓	✓	✓	
184	Weather resistant plastic cabinet with heating (24V)	✓	✓	✓	✓



METAL INDOOR CABINET

High quality

- Strong metal cabinet creates a seamless look
- Ideal for office buildings and public spaces
- Perfectly powder coated
- Logo engraved aluminum hand grip

Dimensions

- 38 x 38 x 20 cm
- Fits defibrillators of all brands

Options

- Acoustic alarm
- Available in green and white



Product no.	Description	Alarm	Fits defibrillators of all brands
250	White metal cabinet		✓
251	White metal cabinet with alarm	✓	✓
252	Green metal cabinet		✓
253	Green metal cabinet with alarm	✓	✓

PERSPEX WALL MOUNT

High quality

- Made of high quality Perspex
- Easy to mount
- Polymer length optimized for both transparency and sturdiness
- Mounting materials included

Dimensions

- 24 x 46 x 15 cm
- Specially designed for HeartSine AED's and other compact defibrillators

Options

- Customized AED stickers



Product no.	Description
37	Perspex wall mount



AED BACKPACK

High quality

- ♥ Designed in collaboration with first-aid services
- ♥ Weather and water resistant
- ♥ Robust character
- ♥ High visibility

Dimensions

- ♥ Compact model: 31 x 40 x 15 cm
- ♥ Large model: 38 x 56 x 19 cm

Options

- ♥ Both models are able to equip defibrillators of all brands
- ♥ Fast response time in case of an emergency
- ♥ Flexible storing solutions



Product no.	Description	Weight	Capacity in litres	Fits defibrillators of all brands	Weather and water resistant
480	Compact AED backpack	1.5 kg	17,5 L	✓	✓
481	Large AED backpack	2.5 kg	35 L	✓	✓



SAFESET

Contains

- Microfiber tissue
- Rescue blanket
- Pocket mask
- Scissors
- Razor
- Gloves



Purpose

- Essential accessories for CPR

Product no.	Description
1	Safe Set

Why choose ARKY?

- ✓ Attractive price-quality ratio
- ✓ High quality materials and components
- ✓ All products are in stock, shipment within 24 hrs
- ✓ Cabinets fit defibrillators of all brands
- ✓ Customization/branding options
- ✓ All products are individually packed in protective plastic and a three-layer carton box

Verlengde Bremenweg 12

9723 JV Groningen

The Netherlands

☎ +31 (0) 50 311 6085

✉ sales@arkycase.com

🌐 www.arkycase.com