

mediTECH

MEDICAL GAS EQUIPMENT

PRODUCT CATALOGUE 2023
END USER PRICE LIST



Reliability and Performance Guaranteed
Manufactured in the United Kingdom

CONTENTS

An Introduction.....	3
Microvents.....	4-8
Regulators.....	9-15
Demand Valves.....	16-17
Flow Selectors.....	18-23
TinyFlow.....	24
Gas Specific Adaptors.....	25 - 26
Hose Assemblies.....	27 - 33
Fittings.....	34 - 40
Complete Ventilator System.....	41
Bedside Oxygen System.....	42
Kits.....	43 - 48
Cylinders.....	49 - 52
Accessories.....	53 - 56
RedVac.....	57 - 58
MedVac.....	59 - 60
Servicing.....	61
Terms and Conditions.....	62-62



AN INTRODUCTION

B.N.O.S. Meditech is a world – class British medical device company, established in Essex in 1989. Meditech designs, manufactures and exports ventilation, resuscitation, therapy + other medical gas products to over one hundred countries worldwide.

We set standards for saving peoples lives. We work closely with professional users within the emergency medical services, hospitals, military medical services, industrial safety + marine to develop market leading medical devices.

Meditech is the UK's leading company in the field of medical gas equipment. Committed to customer service we supply high calibre products and our technical engineering ability permits necessary product modifications to comply with local market requirements throughout the world.

OUR PRODUCTS

B.N.O.S. Meditech Limited manufacturers and distributes products for ventilation, resuscitation, oxygen therapy, pain relief + more. Meditech is an ISO 13485 registered company. All of our products meet the relevant national and international standards and are CE marked in accordance with the Medical Devices Directive in the EU.

In – house manufacture and design give Meditech significant advantage over our competitors. We have extensive CNC machining + injection moulding centres with all other processes of manufacture + assembly performed in – house. Our in – house capability gives Meditech total control of quality, precision, cost + eliminates supply chain.

VISION + VALUES

QUALITY

High quality products + services are at the forefront of everything we do at Meditech. Quality is never compromised.

CUSTOMER SERVICE

We develop long lasting relationships with our customers delivering exceptional pre- and post-market sales. We work together through our shared values, open communication, teamwork + high quality product portfolio.

INNOVATION

We pride ourselves on delivering market leading, innovative products and ensuring we maintain our requirements to meet our customers needs.

SUSTAINABILITY + SOCIAL RESPONSIBILITY

B.N.O.S. Meditech has selected three focus areas for the Social Responsibility policy and sustainability efforts:

- Safety: of our products and in the workplace
- Ethical workplace and business practice
- Usage and recycling of material



MICROVENTS



MICROVENT™ WORLD, ADVANCED, AIRMIX, ADULT AND CHILD, BS SCHRADER INLET FITTING

Part No. 670-0531-00

Price: £1,882.10

Microvent™ resuscitator fitted with:

- Airmix for 50% or 100% Oxygen
- Patient assist valve for breathing synchronisation
- Adjustable selection control for use with adults and children varies the frequency and tidal volume with a 1:21: E ratio
- Automatic (time cycled) or manual operation.

STANDARD option is also available, excluding patient assist functions:

Part No. 670-0539-00

Price: £1,452.00



MICROVENT™ WORLD, ADVANCED, ADULT AND CHILD, BS SCHRADER INLET FITTING

Part No. 670-0506-00

Price: £1,600.23

Microvent™ resuscitator fitted with:

- Patient assist valve for breathing synchronisation
- Adjustable selection control for use with adults and children varies the frequency and tidal volume with a 1:21: E ratio
- Automatic (time cycled) or manual operation.

STANDARD option is also available, excluding patient assist functions:

Part No. 670-0511-00

Price: £1,166.00



MICROVENT™ EUROPEAN, ADVANCED, ADULT AND CHILD, BS SCHRADER INLET FITTING

Part No. 670-0215-00

Price: £1,600.23

Microvent™ resuscitator fitted with:

- Patient assist valve for breathing synchronisation
- Adjustable selection control for use with adults and children varies the frequency and tidal volume with a 1:21: E ratio
- Automatic (time cycled) or manual operation.

STANDARD option is also available, excluding patient assist functions:

Part No. 670-0213-00

Price: £1,166.00

Each Microvent™ is available with the hose fittings shown below. If the hose fitting you require is not shown, please contact us to discuss your requirements.



BS SCHRADER



DISS



DIN



AGA



AFNOR



SISS



CZECH



UNIFOR



MICROVENTS



MICROVENT™ CPR, ADULT AND CHILD, BS SCHRADER INLET FITTING

Part No. 670-0483-00

Price: £1,166.00

Microvent™ resuscitator:

- Set to 10 breaths per minute with a 1 second inspiratory time
- Pneumatic Oxygen powered resuscitator
- Automatic (time cycled) and manual trigger
- Easy to use, lightweight and compact
- Constructed within a tough reinforced plastic casing
- Microvent™ patient valve feature 22/15mm connects to take standard airway devices and face masks
- 2M supply hose with BS Schrader
(other international fittings available)
- Excellent visibility of patient and resuscitator controls



MICROVENT™ INDUSTRIAL, ADVANCED, ADULT ONLY, BS SCHRADER INLET FITTING

Part No. 670-0261-00

Price: £1,237.50

Microvent™ automatic and manual resuscitator fitted with:

- Patient assist valve for breathing synchronisation
- Designed for easy use in industrial settings
- Manual setting ideal for 30:2 ventilation in CPR

STANDARD option is also available, excluding patient assist functions:

Part No. 670-0259-00

Price: £808.50



MICROVENT™ UTILITY, ADVANCED, ADULT, BS SCHRADER INLET FITTING

Part No. 670-0339-00

Price: £1,144.55

The "switch and use" automatic resuscitator:

- Patient assist valve for breathing synchronisation
- Designed for easy use in industrial settings
- Manual setting ideal for 30:2 ventilation in CPR

STANDARD option is also available, excluding patient assist functions:

Part No 670-0583-00

Price: £715.55

Each Microvent™ is available with the hose fittings shown below. If the hose fitting you require is not shown, please contact us to discuss your requirements.



BS SCHRADER



DISS



DIN



AGA



AFNOR



SISS



CZECH



UNIFOR



MICROVENTS



MICROVENT™ RESPONDER RESUSCITATOR, BS SCHRADER INLET FITTING

Part No. 670-0312-00

Price: £379.50

Microvent™ Responder resuscitator fitted with:

- Manual control
- Flow rates of 20lpm or 40lpm are selected
- For Adult & Child use
- Trigger operated to start ventilation
- Controlled flow ventilation (with pressure relief) and no Oxygen wastage.



FABRIC BAG AND ACCESSORIES FOR PROFESSIONAL MODELS FOR ADULTS AND CHILDREN

Part No. 131-0069-20

Price: £141.90

Wipeable carrying bag with shoulder strap, carrying handles, internal pockets, cylinder restraint straps, etc. Supplied with:

- 3 X Resuscitation masks (adult, small adult, child)
- 2 X High concentration therapy masks (adult, child)
- 4 X Oropharyngeal airways (size 0, 1,2,3)
- 1 X ET tube connector

****Available at this price only when purchased with a Microvent™ resuscitator ****

FABRIC BAG AND ACCESSORIES FOR INDUSTRIAL AND UTILITY MODEL MICROVENT™ FOR ADULTS

Part No. 131-0079-20

Price: £123.75

Wipeable carrying bag with shoulder strap, carrying handles, internal pockets, cylinder restraint straps, etc. Supplied with:

- 1 X Resuscitation mask (adult)
- 1 X High concentration therapy mask (adult)
- 2 X Oropharyngeal airways
- 1 X ET tube connector

****Available at this price only when purchased with a Microvent™ resuscitator ***



MICROVENT™ CIRCUIT INCLUDING PEEP

Part No.

Price

673-1040-00

£38.50

- 22MM Flextube transport breathing system with non-re breathing valve + adjustable PEEP valve. 0-20CM H₂O (Fixed PEEP available upon request)
- Includes Microvent™ protection filter, HME filter, adult + child resuscitation masks.

MICROVENT™ CIRCUIT EXCLUDING PEEP

Part No.

Price

673-1040-01

£35.20

- 22MM Flextube transport breathing system with non-re breathing valve
- Includes Microvent™ protection filter, HME filter, adult + child resuscitation masks.



Child

Small Adult

Adult

SINGLE USE RESUSCITATION FACE MASKS

Size	Part No.	Price
Child	673-0558-00	£16.50
Small Adult	673-0559-00	£16.50
Adult	673-0582-00	£16.50

Single use face masks. Anatomical shape for a tailored fit.


- Soft seal for comfort.
- Transparent shell for excellent visibility.
- Manufactured from polypropylene and thermoplastic elastomer, contains no PVC and is phthalate-free.



HOOK RING - For microVENT

Part No.	Price
673-0444-00	£3.08

The hook ring attaches the head harness to the silicone resuscitation mask. For use with head harness 673-0445-00 and resuscitation face mask 673-0452-00.



HEAD HARNESS - For microVENT

Part No.	Price	Part No.	Price
673-0445-00 (Large adult)	£9.46	673-0445-01 (Small adult)	£9.46

The head harness holds the resuscitation face mask to the patient's mouth and nose, such as in an emergency rescue. For use with hook ring 673-0444-00 and resuscitation face mask 673-0452-00.



I-GEL® SUPRAGLOTTIC AIRWAYS

Size	Part No.	Price
5	673-0000-00	£11.90
4	673-0000-01	£11.90
3	673-0000-02	£11.90
2	673-0000-03	£11.90
1	673-0000-04	£11.90

- Supraglottic airway device for use with Microvent™ Resuscitator
- Designed to create a non-inflatable, anatomical seal of the pharyngeal, laryngeal and perilaryngeal structures. Which avoiding compression trauma
- Ideal for use with patient weights between 2-90KG



CLEAR-THERM 3 FILTER / HME

Part No.	Price
673-1041-00	£4.13

- 22F / 15M – 22M – 15F connectors
- Filtration efficiency is more than 99.9%
- Moisture loss 7.8mg H₂O / l



CG MIDI BREATHING FILTER

Part No.	Price
673-1042-00	£4.13

- 22F – 22M / 15F connectors
- Filtration efficiency is more than 99.9%



CATHETER MOUNT

Part No.	Price
673-0310-00	£3.74

ET tube or resuscitation mask extension mount. Single patient use. Colour and detail may vary.



PRESSURE INDICATOR ON 'T' 720 SWIVEL -For microVENT

Part No.	Price
674-0018-00	£4.18

Pressure indicator with a 22mm male/22mm female fitting. Swivels 360 by 360. Low cost for single use.



EXPIRATION DIVERTER FOR MICROVENT™/ ENTONOX DEMAND VALVE SCAVAGING PORT

Part No.	Price
674-0008-00	£14.19

PEEP valve 674-0007-00 attaches to Microvent™ resuscitator via the expiration diverter 674-0008-00.



PEEP VALVE FOR MICROVENT™

Part No.	Price
674-0007-00	£20.13

PEEP valve 674-0007-00 attaches to Microvent™ resuscitator via the expiration diverter 674-0008-00.

OXYGEN REGULATOR, PIN INDEX, THERAPY ONLY.



Part No.	Price
55B 677-0000-00	£163.35
55SL 777-0000-00	£163.35

- Multiple flow therapy outlet set by selector dial – ½, 1,2,3,4,6,8,10 and 15 litres per minute.
- 0-25 & 0-70lpm flow rates also available
- 55B approx. 634g weight
- 55SL approx. 385g weight

OXYGEN REGULATOR, PIN INDEX, 2 x BS SCHRADER THERAPY.



Part No.	Price:
55B 677-0001-20	£254.10
55SL 777-0001-20	£254.10

- Multiple flow therapy outlet set by selector dial – ½, 1,2,3,4,6,8,10 and 15 litres per minute.
- 0-25 & 0-70lpm flow rates also available
- 55B approx. 847g weight
- 55SL approx. 538g weight

OXYGEN REGULATOR, PIN INDEX, 2 x BS SCHRADER.



Part No.	Price
55B 678-0001-20	£206.25
55SL 778-0001-20	£206.25

- Non-therapy Regulator
- 55B approx. 776g weight
- 55SL approx. 469g weight

OXYGEN REGULATOR, PIN INDEX, 1 x BS SCHRADER THERAPY.



Part No.	Price
55B 677-0001-10	£211.75
55SL 777-0001-10	£211.75

- Multiple flow therapy outlet set by selector dial – ½, 1,2,3,4,6,8,10 and 15 litres per minute.
- 0-25 & 0-70lpm flow rates also available
- 55B approx. 795g weight
- 55SL approx. 485g weight

OXYGEN REGULATOR, PIN INDEX, 1 x BS SCHRADER.



Part No.	Price
55B 678-0001-10	£163.35
55SL 778-0001-10	£163.35

- Non-therapy Regulator
- 55B approx. 724g weight
- 55SL approx. 416g weight

- 55B = Brass body
- 55SL = Super light body

Available Gases:

- Oxygen (O₂)
- Medical Air (MA4)
- Surgical Air (MA7)
- Nitrous Oxide (N₂O)
- Entonox (N₂O/O₂)
- Carbon Dioxide (CO₂)

MEDICAL AIR REGULATOR, PIN INDEX, THERAPY ONLY.



Part No.	Price
55B 691-0900-00	£163.35
55SL 791-0900-00	£163.35

- Multiple flow therapy outlet set by selector dial – ½, 1,2,3,4,6,8,10 and 15 litres per minute.
- 0-25 & 0-70lpm flow rates also available
- 55B approx. 634g weight
- 55SL approx. 385g weight

MEDICAL AIR REGULATOR, PIN INDEX, 2 X BS SCHRADER, THERAPY.



Part No.	Price
55B 691-0908-20	£254.65
55SL 791-0908-20	£254.65

- Multiple flow therapy outlet set by selector dial – ½, 1,2,3,4,6,8,10 and 15 litres per minute.
- 0-25 & 0-70lpm flow rates also available
- 55B approx. 847g weight
- 55SL approx. 538g weight

MEDICAL AIR REGULATOR, PIN INDEX, 2 X BS SCHRADER.



Part No.	Price
55B 678-0001-20	£206.25
55SL 778-0001-20	£206.25

- Non-therapy Regulator
- 55B approx. 776g weight
- 55SL approx. 469g weight

MEDICAL AIR REGULATOR, PIN INDEX, 1 X BS SCHRADER, THERAPY.



Part No.	Price
55B 691-0908-10	£211.75
55SL 791-0908-10	£211.75

- Multiple flow therapy outlet set by selector dial – ½, 1,2,3,4,6,8,10 and 15 litres per minute.
- 0-25 & 0-70lpm flow rates also available
- 55B approx. 795g weight
- 55SL approx. 485g weight

MEDICAL AIR REGULATOR, PIN INDEX, 1 X BS SCHRADER.



Part No.	Price
55B 692-0908-10	£163.35
55SL 792-0908-10	£163.35

- Non-therapy Regulator
- 55B approx. 724g weight
- 55SL approx. 416g weight

- 55B = Brass body
- 55SL = Super light body

Available Gases:

- Oxygen (O₂)
- Medical Air (MA4)
- Surgical Air (MA7)
- Nitrous Oxide (N₂O)
- Entonox (N₂O/O₂)
- Carbon Dioxide (CO₂)

OXYGEN REGULATOR, BULLNOSE INLET, THERAPY ONLY.



Part No.	Price
55B 679-0100-00	£163.35
55SL 779-0100-00	£163.35

- Multiple flow therapy outlet set by selector dial – ½, 1,2,3,4,6,8,10 and 15 litres per minute.
- 0-25 & 0-70lpm flow rates also available
- 55B approx. 668g weight
- 55SL approx. 485g weight

OXYGEN REGULATOR, BULLNOSE INLET, 2 X BS SCHRADER, THERAPY.



Part No.	Price
55B 679-0101-20	£254.10
55SL 779-0101-20	£254.10

- Multiple flow therapy outlet set by selector dial – ½, 1,2,3,4,6,8,10 and 15 litres per minute.
- 0-25 & 0-70lpm flow rates also available
- 55B approx. 772g weight
- 55SL approx. 588g weight

OXYGEN REGULATOR, BULLNOSE INLET, 2 X BS SCHRADER.



Part No.	Price
55B 680-0101-20	£206.25
55SL 780-0101-20	£206.25

- Non-therapy Regulator
- 55B approx. 772g weight
- 55SL approx. 588g weight

OXYGEN REGULATOR, BULLNOSE INLET, 1 X BS SCHRADER, THERAPY.



Part No.	Price
55B 679-0101-10	£211.75
55SL 779-0101-10	£211.75

- Multiple flow therapy outlet set by selector dial – ½, 1,2,3,4,6,8,10 and 15 litres per minute.
- 0-25 & 0-70lpm flow rates also available
- 55B approx. 720g weight
- 55SL approx. 537g weight

OXYGEN REGULATOR, BULLNOSE INLET, 1 X BS SCHRADER.



Part No.	Price
55B 680-0101-10	£163.35
55SL 780-0101-10	£163.35

- Non-therapy Regulator
- 55B approx. 645g weight
- 55SL approx. 463g weight

TOP ENTRY BULLNOSE INLET REGULATOR, THERAPY ONLY



Part No.	Price
55B 673-0646-00	£206.25
55SL 773-0646-00	£206.25

- Top Entry Oxygen Regulator
- 55B approx. 668g weight
- 55SL approx. 485g weight

MEDICAL AIR REGULATOR, BULLNOSE INLET, THERAPY ONLY.



Part No.	Price
55B 693-0100-00	£163.85
55SL 793-0100-00	£163.85

- Multiple flow therapy outlet set by selector dial – ½, 1,2,3,4,6,8,10 and 15 litres per minute.
- 0-25 & 0-70lpm flow rates also available
- 55B approx. 668g weight
- 55SL approx. 485g weight

MEDICAL AIR REGULATOR, BULLNOSE INLET, 2 X BS SCHRADER, THERAPY.



Part No.	Price
55B 693-0108-20	£254.10
55SL 793-0108-20	£254.10

- Multiple flow therapy outlet set by selector dial – ½, 1,2,3,4,6,8,10 and 15 litres per minute.
- 0-25 & 0-70lpm flow rates also available
- 55B approx. 772g weight
- 55SL approx. 588g weight

MEDICAL AIR REGULATOR, BULLNOSE INLET, 2 X BS SCHRADER.



Part No.	Price
55B 694-0108-20	£206.75
55SL 794-0108-20	£206.75

- Non-therapy Regulator
- 55B approx. 772g weight
- 55SL approx. 588g weight

MEDICAL AIR REGULATOR, BULLNOSE INLET, 1 X BS SCHRADER, THERAPY.



Part No.	Price
55B 693-0108-10	£211.75
55SL 793-0108-10	£211.75

- Multiple flow therapy outlet set by selector dial – ½, 1,2,3,4,6,8,10 and 15 litres per minute.
- 0-25 & 0-70lpm flow rates also available
- 55B approx. 720g weight
- 55SL approx. 537g weight

MEDICAL AIR REGULATOR, BULLNOSE INLET, 1 X BS SCHRADER.



Part No.	Price
55B 694-0108-10	£163.85
55SL 794-0108-10	£163.85

- Non-therapy Regulator
- 55B approx. 645g weight
- 55SL approx. 463g weight

- 55B = Brass body
- 55SL = Super light body

Available Gases:

- Oxygen (O₂)
- Medical Air (MA4)
- Surgical Air (MA7)
- Nitrous Oxide (N₂O)
- Entonox (N₂O/O₂)
- Carbon Dioxide (CO₂)

NON-MAGNETIC PIN INDEX MEDICAL GAS REGULATOR, THERAPY ONLY



Part No.	Price
720-0000-00 (0-15lpm)	£248.05
720-0000-01 (0-25lpm)	£248.05

- Suitable for MRI scans
- Non-magnetic material
- No calibration
- Available with a Pin index, Bullnose & CGA-540 inlet

NON-MAGNETIC PIN INDEX MEDICAL GAS REGULATOR, 1 x BS SCHRADER, THERAPY



Part No.	Price
720-0001-11 (0-15lpm)	£248.05
720-0001-01 (0-25lpm)	£248.05

- Suitable for MRI scans
- Non-magnetic material
- No calibration
- Available with a Pin index, Bullnose & CGA-540 inlet

NITROUS OXIDE REGULATOR WITH BULLNOSE B5341 NO.13 INLET, 1 X BS SCHRADER OUTLET



Part No.	Price
673-0114-10	£203.50

- 5-year warranty
- 5-year service interval
- Available with any international HP outlet fitting
- Available with a Pin index, Bullnose / B5341 NO.13 inlet

NITROUS OXIDE REGULATOR, WITH PIN INDEX INLET, 1 X BS NITROUS OXIDE BS SCHRADER OUTLET



Part No.	Price
55B 699-0014-10	£163.85
55SL 799-0014-10	£163.85

- 5-year warranty
- 5-year service interval
- Available with any international HP outlet fitting
- Available with a Pin index, Bullnose / B5341 NO.13 inlet

ENTONOX® REGULATOR WITH PIN INDEX INLET, 1 X ENTONOX® BS SCHRADER OULET



Part No.	Price
55B 673-0546-00	£163.85
55SL 673-0546-01	£163.85

- 5-year warranty
- 5-year service interval
- Available with any international HP outlet fitting
- Available with a Pin index, Bullnose / B5341 NO.13 inlet

CO2 REGULATOR WITH PIN INDEX INLET, 1 X CO2 BS SCHRADER OULET



Part No.	Price
55B 673-0433-00	£163.85
55SL 673-0419-00	£163.85

- 5-year warranty
- 5-year service interval
- Available with any international HP outlet fitting
- Available with a Pin index, Bullnose / B5341 NO.13 inlet

Create a regulator to suit your specific requirements by following the steps below to generate your regulator part number.

STEP 1 Choose the Regulator **BODY TYPE** that you require. Decide if you want Brass or Aluminium.

STEP 2 Choose the **TYPE** of Oxygen Pressure **INLET** that you require.

STEP 3 Choose the **TYPE** of Oxygen Pressure **OUTLET** that you require.

STEP 4 Choose **HOW MANY** Oxygen Pressure Outlets that you need for your Regulator.

STEP 5 Decide if you want your Regulator to have **OXYGEN THERAPY**, and if so, choose what **TYPE**.

Basic Regulator Type:

BRASS - 55B

- 677 – Pin Index regulator with therapy
- 678 – Pin Index regulator with no therapy
- 679 – 55B International (nut fitting) therapy regulator
- 680 – 55B International (nut fitting) non-therapy regulator

ALUMINIUM – 55SL

- 777 – Pin index regulator with therapy
- 778 – Pin index without therapy
- 779 – 55SL international (nut fitting) therapy regulator
- 780 – 55SL international (nut

4 BAR MEDICAL GAS REGULATOR PART NUMBERING SYSTEM

Type of Outlet:

- 00 = (therapy only)
- 01 = BS5682 BS Schrader (UK)
- 02 = DISS 1240 (USA)
- 03 = AGA SS875 24 30:2004
- 04 = DIN 13260-2:2004
- 05 = AFNOR NF S90-116:1988 (French)
- 07 = SIS (sleeve index system) AS2902:2016 (Australian)
- 08 = MA4
- 19 = HOSE BARB 1/4"

Number of Outlet Fittings:

- 0 = None (Therapy ONLY)
 - 1 = One (to Front)
 - 2 = Two
 - 3 = One (to Rear)
- NOTE: Front Outlet is adjacent to Gauge



To create a Pin-Index Articulating Yoke Regulator, start the Product Code with

- 679-00 - 55B**
- 779-00 - 55SL**

Type of inlet:

- 00 = Pin index BS EN ISO407:2021 (CGA 870, USA)
- 01 = Bullnose (UK)
- 02 = Bullnose (Benelux)
- 03 = W21.8 1/14th
- 04 = DINN 477 (G 3/4")
- 05 = CGA540 (USA)

Maximum Therapy Flow Rate:

- 0 = Therapy flows 0-15lpm (fir tree connector)
- 0 = No therapy (use 0 also)
- 1 = Therapy flows 0-25lpm (fir tree connector)
- 2 = Therapy flows 0-15lpm (Threaded connector 9/16" – 18 for humidifier)
- 3 = Therapy flows 0-25lpm (Threaded connector 9/16" – 18 for humidifier)
- 4 = 1/4 "HOSE BARB therapy flows 0-15lpm
- 5 = 1/4 "HOSE BARB therapy flows 0-25lpm
- 6 = therapy flow 0-70lpm (threaded connector 9/16" – 18 for humidifier)
- 7 = 0 – 70 lpm (fir tree)
- 8 = 3/8 BSP 0 – 15 lpm
- 9 = 3/8 BSP 0-25 lpm



Articulating Yoke

INLET TYPES



PIN INDEX



BULLNOSE



BENELUX BULLNOSE



W21.8 1/2"



DIN 477 (G 3/4")



CGA540



Pin Index Yoke
Articulating



OUTLET TYPES



BS SCHRADER



DISS



DIN



AGA



AFNOR



SIS



CZECH



UNIFOR



THERAPY TYPES



FIR TREE



3/8" THREADED
THERAPY CONNECTOR



3/8" BSP THREADED
THERAPY CONNECTOR



1/4" BARB

For further information or requirements regarding our International Regulators, please contact a member of the sales team.

DEMAND VALVES

ENTONOX® DEMAND VALVE



Part No.	Price
900-0021-00	£260.15 (2M HOSE)
900-0031-00	£272.25 (3M HOSE)
900-0041-00	£284.35 (4M HOSE)

- Available in 2M,3M & 4M hose with BS Schrader probe
- Full length handle
- 5 Year Warranty
- 5 Year service
- 10 year intended life-span

Our analgesic demand valve is designed for use with 50/50 N₂O/O₂ (Nitronox/Entonox®) to meet BS4272-2. It is extremely strong with practical features to avoid cross contamination which provides a good patient experience. There is low patient effort required. Our demand valves are designed to use with generic filter systems. Each demand valve can be supplied with a dedicated pressure regulator or with any international hose fitting to use with your own cylinder or wall supply.



ENTONOX® REGULATOR PIN INDEX 1 x BS SCHRADER

Part No.	Price
673-0546-00	£163.35

- Entonox® regulator with pin index inlet
- 1 x BS Schrader outlet

BACTERIAL FILTER/MOUTH PIECE

Part No.	Price
673-9021-00	£8.50

Single use bacterial mouth pieces used to avoid cross contamination. To be used with our Entonox® and Oxygen demand valves.

EXPIRATION AGSS DIVERTER

Part No.	Price
674-0008-00	£12.90

Expiration AGSS Diverter attaches to Demand Valve for scavaging exhaled gas.



WRIST BAND

Part No.	Price
674-9000-00	£27.50

Each Demand Valve is available with the hose fittings shown below. If the hose fitting you require is not shown, please contact us to discuss your requirements.



BS SCHRADER



DISS



DIN



AGA



AFNOR



SISS



CZECH



UNIFOR



DEMAND VALVES

OXYGEN DEMAND VALVE

Part No.	Price
910-0021-00	£260.15 (2M HOSE)
910-0031-00	£272.25 (3M HOSE)
910-0041-00	£284.35 (4M HOSE)

- Available in 2M,3M & 4M hose with BS Schrader probe
- Full length handle
- 5 Year Warranty
- 5 Year service
- 10 year intended life-span

Our Oxygen demand valve is designed for use with 100% O₂ to meet BS4272-2. It is extremely strong with practical features to avoid cross contamination which provides a good patient experience. There is low patient effort required. Our demand valves are designed to use with generic filter systems. Each demand valve can be supplied with a dedicated pressure regulator or with any international hose fitting to use with your own cylinder or wall supply.

OXYGEN DEMAND VALVE - NO HANDLE

Part No.	Price
910-0027-00	£242.00

Our Oxygen Demand Valve with no handle has been designed to work with BVM resuscitators to eliminate any wasted Oxygen.



OXYGEN REGULATOR – PIN INDEX 1 x BS SCHRADER

Part No.	Price
678-0001-10	£163.35

DEMAND VALVE TO BVM ADAPTOR

AMBU Part No.	Price	LAERDAL Part No.	Price
673-0636-00	£21.50	673-0673-00	£21.50

Each Demand Valve is available with the hose fittings shown below. If the hose fitting you require is not shown, please contact us to discuss your requirements.



BS SCHRADER



DISS



DIN



AGA



AFNOR



SISS



CZECH



UNIFOR



FLOW SELECTORS



55SL FLOW SELECTOR, BS SCHRADER PROBE

lpm	Part No.	Price
0-15	673-0506-00	£120.45
0-25	673-0506-01	£120.45
0-70	673-0070-00	£132.55

- Flow selector with multiple flow therapy outlet set by selector dial – ½, 1,2,3,4,6,8,10,15,25 and 70 litres per minute.
- Designed for use with a 4 bar Oxygen supply
- Available with left angle probe.
- Available with any international probe type.



55SL FLOW SELECTOR WITH EXTENDED BS SCHRADER PROBE

lpm	Part No.	Price
0-15	673-0585-00	£127.60
0-25	673-0585-01	£127.60
0-70	673-0071-00	£132.55

- Flow selector with multiple flow therapy outlet set by selector dial – ½, 1,2,3,4,6,8,10,15,25 and 70 litres per minute.
- Designed for use with a 4 bar Oxygen supply
- Available with left angle probe.
- Available with any international probe type.



TWIN FLOW SELECTOR WITH 2M SUPPLY HOSE FITTED WITH BS SCHRADER PROBE PLUS 2X HUMIDIFIERS

lpm	Part No.	Price
0-15	710-0101-22	£290.40*
0-25	710-0101-23	£290.40*
0-70	710-0101-24	£302.50*

- Twin outlet fittings are 9/16" thread to take humidifiers
 - Converters are also supplied to take Oxygen therapy hose
 - The selector has a wall mounting plate allowing fixing to a wall or bulkhead by screws or bolts.
 - Available with any international probe type.
- *Humidifiers to be purchased separately.*

TWIN FLOW SELECTOR RAIL MOUNT WITH 2M SUPPLY HOSE FITTED BS SCHRADER PROBE PLUS 2X HUMIDIFIERS

lpm	Part No.	Price
0-15	710-0102-22	£326.70*
0-25	710-0102-23	£326.70*
0-70	710-0102-24	£338.80*

- Twin outlet fittings are 9/16" thread to take humidifiers
 - Converters are also supplied to take Oxygen therapy hose
 - The selector has a rail mounting plate allowing fixing to a medical rail.
 - Available with any international probe type.
- *Humidifiers to be purchased separately.*

Available gas Types:

- Oxygen (O₂)
- Medical Air (MA4)
- Entonox® (N₂O/O₂)

THERAPY TYPES



FIR TREE



3/16" THREADED
THERAPY
CONNECTOR



3/8" BSP THREADED
THERAPY
CONNECTOR



1/4" BARB

FLOW SELECTORS

55SL 4 WAY FLOW SELECTOR AND FRAME, BS SCHRADER PROBE ON 2M SUPPLY HOSE



lpm	Part No.	Price
0-15	673-0534-00	£1,028.50
0-25	673-0534-01	£1,028.50
0-70	673-0534-07	£1,040.60

- 55SL 4-way flow selector with frame
- Flow selector with multiple flow therapy outlet set by selector dial – ½, 1,2,3,4,6,8,10,15,25 and 70 litres per minute.
- Designed for use with a 4 bar Oxygen supply
- Approximate dimensions 460 X 210 X 210mm
- Also available with different international hose fittings.

55SL 4-WAY FLOW SELECTOR ON STAND-MOUNT, BS SCHRADER PROBE, 2M SUPPLY HOSE



lpm	Part No.	Price
0-15	673-0834-00	£665.50
0-25	673-0835-00	£665.50
0-70	673-0870-00	£677.60

- Flow selector with multiple flow therapy outlet set by selector dial – ½, 1,2,3,4,6,8,10,15,25 and 70 litres per minute.
- Designed for use with a 4 bar Oxygen supply
- Also available with left angle probe. Please contact us for details.
- Also available with different international hose fittings.

55SL 4-WAY FLOW SELECTOR ON STAND-MOUNT, BS CHECK UNIT



lpm	Part No.	Price
0-15	673-0934-00	£665.50
0-25	673-0935-00	£665.50
0-70	673-0970-00	£677.60

- Flow selector with multiple flow therapy outlet set by selector dial – ½, 1,2,3,4,6,8,10,15,25 and 70 litres per minute.
- Designed for use with a 4 bar Oxygen supply
- Also available with different international hose fittings.

Available gas Types:

- Oxygen (O₂)
- Medical Air (MA4)
- Entonox® (N₂O/O₂)

THERAPY TYPES



FIR TREE



3/8" THREADED
THERAPY
CONNECTOR



3/8" BSP THREADED
THERAPY
CONNECTOR



1/4" BARB

FLOW SELECTORS



Twin Flow Selector, with direct BS Probe

lpm	Part No.	Price
0-15	710-0106-22	£200.75
0-25	710-0106-23	£200.75
0-70	710-0106-24	£212.85

- Twin Flow Selector with direct BS Probe
- Flow selector with multiple flow therapy outlet set by selector dial – ½, 1,2,3,4,6,8,10,15,25 and 70 litres per minute.
- Designed for use with a 4 bar Oxygen supply
- Also available with different international probe types



Rail Mounted Single Flow Selector, with 2M Hose Assembly, BS Probe

lpm	Part No.	Price
0-15	673-2003-00	£179.30
0-25	673-2004-00	£179.30
0-70	673-2001-00	£190.30

- Single Flow Selector, Rail Mounted, with 2M supply hose
- Flow selector with multiple flow therapy outlet set by selector dial – ½, 1,2,3,4,6,8,10,15,25 and 70 litres per minute.
- Designed for use with a 4 bar Oxygen supply
- Also available with different international probe types



Rail Mounted Twin Flow Selector, with 2M Hose Assembly, BS Probe

lpm	Part No.	Price
0-15	710-0107-22	£281.49
0-25	710-0107-23	£281.49
0-70	710-0107-24	£292.49

- Twin Flow Selector, Rail Mounted, with 2M supply hose
- Flow selector with multiple flow therapy outlet set by selector dial – ½, 1,2,3,4,6,8,10,15,25 and 70 litres per minute.
- Designed for use with a 4 bar Oxygen supply
- Also available with different international probe types



Twin Flow Selector, Twin Adaptor with BS Probe

lpm	Part No.	Price
0-15	710-0108-22	£270.49
0-25	710-0108-23	£270.49
0-70	710-0108-24	£281.49

- Twin Flow Selector, Twin Adaptor with direct BS Probe
- Flow selector with multiple flow therapy outlet set by selector dial – ½, 1,2,3,4,6,8,10,15,25 and 70 litres per minute.
- Designed for use with a 4 bar Oxygen supply
- Also available with different international probe types

Available gas Types:

- Oxygen (O₂)
- Medical Air (MA4)
- Entonox (N₂O/O₂)

THERAPY TYPES



FIR TREE



3/8" THREADED THERAPY CONNECTOR



3/8" BSP THREADED THERAPY CONNECTOR



1/4" BARB

FLOW SELECTORS



Twin Flow Selector, Straight Bar – 100mm – BS Probe

lpm	Part No.	Price
0-15	710-0109-22	£270.49
0-25	710-0109-23	£270.49
0-70	710-0109-24	£281.49

- Twin Flow Selector mounted on a 100mm straight bar with BS Probe
- Flow selector with multiple flow therapy outlet set by selector dial – ½, 1,2,3,4,6,8,10,15,25 and 70 litres per minute.
- Designed for use with a 4 bar Oxygen supply
- Also available with different international probe types.



Twin Flow Selector, Rail Mounted, Straight Bar – 100mm – 2M Supply Hose, BS Probe

lpm	Part No.	Price
0-15	710-0110-22	£281.49
0-25	710-0110-23	£281.49
0-70	710-0110-24	£292.49

- Twin Flow Selector on straight bar on a rail mount – 2M supply hose, BS Probe
- Flow selector with multiple flow therapy outlet set by selector dial – ½, 1,2,3,4,6,8,10,15,25 and 70 litres per minute.
- Designed for use with a 4 bar Oxygen supply
- Also available with different international probe types.



Straight Bar, Flow Selector + BS Schrader Outlet – Direct BS Probe – 100MM

lpm	Part No.	Price
0-15	673-3000-00	£259.49
0-25	673-3001-00	£259.49
0-70	673-3002-00	£270.49

- Straight bar with 1x Flow Selector, 1x BS Schrader Outlet – direct BS Probe – 100MM
- Flow Selector with multiple flow therapy outlet set by selector dial – ½, 1,2,3,4,6,8,10,15,25 and 70 litres per minute.
- Designed for use with a 4 bar Oxygen supply
- Also available with different international probe types.



Straight Bar, Rail Mounted, Flow Selector + BS Schrader Outlet – 2m Supply Hose, BS Probe – 100MM

lpm	Part No.	Price
0-15	673-3003-00	£270.49
0-25	673-3004-00	£270.49
0-70	673-3005-00	£281.49

- Straight bar, rail mounted with 1x Flow Selector, 1x BS Schrader Outlet – 2M supply hose – 100MM
- Flow selector with multiple flow therapy outlet set by selector dial – ½, 1,2,3,4,6,8,10,15,25 and 70 litres per minute.
- Designed for use with a 4 bar Oxygen supply
- Also available with different international probe types.

Available gas Types:

- Oxygen (O₂)
- Medical Air (MA4)
- Entonox (N₂O/O₂)

THERAPY TYPES



FIR TREE



3/8" THREADED THERAPY CONNECTOR



3/8" BSP THREADED THERAPY CONNECTOR



1/4" BARB



Bobbin Style Flow Selector-BS Probe

Gas	lpm	Part No.	Price
O ₂	0-15	812-0000-00	£42.35
O ₂	0-25	812-0001-00	£42.35
O ₂	0-70	812-0002-00	£75.35

- Designed for use with a 4 bar Oxygen supply
- Also available with different international probe types



Bobbin Style Flow Selector- BS Probe

Gas	lpm	Part No.	Price
MA4	0-15	813-0000-00	£42.35
MA4	0-25	813-0001-00	£42.35
MA4	0-70	813-0002-00	£75.35

- Designed for use with a 4 bar Oxygen supply
- Also available with different international probe types



Dual Flow Bobbin Style Flow Selector - BS Probe

Gas	lpm	Part No.	Price
O ₂	0-15	814-0000-00	£107.50
O ₂	0-25	814-0001-00	£107.50
O ₂	0-70	814-0002-00	£163.35

Gas	lpm	Part No.	Price
MA4	0-15	815-0000-00	£107.50
MA4	0-25	815-0001-00	£107.50
MA4	0-70	815-0002-00	£163.35

- Designed for use with a 4 bar Oxygen supply
- Also available with different international probe types

PROBE TYPE



BS PROBE



DISS



DIN



AGA



AFNOR



SIS



CZECH



UNI



OHMEDA



OUTLET TYPES



PURITAN



CHEMETRON



3/8" BSP Nut +
Liner



1/4" THREADED
THERAPY
CONNECTOR



3/8" BSP THREADED
THERAPY
CONNECTOR

Each Bobbin Style Flow Selector is available with the probe types and outlet types shown. If either the probe type or outlet type you require is not shown, please contact us to discuss your requirements.

Create a Bobbin Style Flow Selector to suit your specific requirements by following the steps below to generate your BOBBIN STLYE FLOW SELECTOR PART NUMBER.


STEP 1 Choose the **GAS TYPE** that you require

STEP 2 Choose the **PROBE TYPE** you require

STEP 3 Choose the **FLOW RATE** you require

STEP 4 Choose the **OUTLET TYPE** you require

BOBBIN STYLE PART NUMBERING SYSTEM





TINYFLOW NEONATAL FLOW SELECTOR, BS SCHRADER PROBE

LPM	Part No.	Price	Gas
0-1	673-0503-00	£195.25	O ₂
0-1	673-0503-01	£195.25	MA4

- TinyFlow with multiple micro-flow rates set by selector dial.
- 0.01, 0.02, 0.03, 0.04, 0.05, 0.06, 0.07, 0.08, 0.09, 0.10 & 1 LPM
- Designed for use with a 4 bar Oxygen supply
- For use with neonatals.
- Available with any international probe type.



TINYFLOW PEDIATRIC FLOW SELECTOR, BS SCHRADER PROBE

lpm	Part No.	Price	Gas
0-3	673-0502-00	£195.25	O ₂
0-3	673-1502-00	£195.25	MA4

- TinyFlow with multiple micro-flow rates set by selector dial.
- 0.02, 0.03, 0.05, 0.08, 0.12, 0.20, 0.30, 0.50, 0.75, 1 & 3 LPM
- Designed for use with a 4 bar Oxygen supply
- For use with pediatric.
- Available with any international probe type.



LPM	Part No.	Price
0-1	776-0000-00	£220.00
0-3	776-0000-01	£220.00

MEDICAL OXYGEN TINYFLOW REGULATOR, PIN INDEX, THERAPY ONLY.

- Oxygen TinyFlow Regulator with Pin Index inlet
- TinyFlow with multiple micro-flow rates set by selector dial.
- 0-1 LPM: 0.01, 0.02, 0.03, 0.04, 0.05, 0.06, 0.07, 0.08, 0.09, 0.10 & 1 litres per minute.
- 0-3LPM: 0.02, 0.03, 0.05, 0.08, 0.12, 0.20, 0.30, 0.50, 0.75, 1 & 3 litres per minute.



LPM	Part No.	Price
0-1	776-0100-00	£220.00
0-3	776-0100-01	£220.00

MEDICAL OXYGEN TINYFLOW REGULATOR, BULLNOSE INLET, THERAPY ONLY.

- Oxygen TinyFlow Regulator with Bullnose inlet
- TinyFlow with multiple micro-flow rates set by selector dial.
- 0-1 LPM: 0.01, 0.02, 0.03, 0.04, 0.05, 0.06, 0.07, 0.08, 0.09, 0.10 & 1 litres per minute.
- 0-3LPM: 0.02, 0.03, 0.05, 0.08, 0.12, 0.20, 0.30, 0.50, 0.75, 1 & 3 litres per minute.

GAS SPECIFIC ADAPTORS



Straight Bar Twin Adaptor – BS Schrader – BS Probe

Gas	Part No.	Price
O ₂	673-9025-00	£120.73
MA4	673-9025- 01	£120.73
MA7	673-9025-02	£120.73
O ₂ /N ₂ O	673-9025-03	£120.73
N ₂ O	673-9025-04	£120.73
Vacuum	673-9025-05	£120.73
CO ₂	673-9025-06	£120.73

- Straight Bar Twin Adaptor, 2x BS Schrader Outlets + BS Probe
- 100mm Straight Bar (50mm available upon request)
- Designed for use with a 4 bar Oxygen supply
- Also available with different international outlet + probe types.



Straight Bar Twin Adaptor – BS Schrader – 3/8” Nut + Nipple

Gas	Part No.	Price
O ₂	673-9035-00	£109.73
MA4	673-9035- 01	£109.73
MA7	673-9035-02	£109.73
O ₂ /N ₂ O	673-9035-03	£109.73
N ₂ O	673-9035-04	£109.73
Vacuum	673-9035-05	£109.73
CO ₂	673-9035-06	£109.73

- Straight Bar Twin Adaptor, 2x BS Schrader Outlets + 3/8” BSP Nut + Nipple
- 100mm Straight Bar (50mm available upon request)
- Designed for use with a 4 bar Oxygen supply
- Also available with different international outlet + probe types.

Each Gas Specific Adapter is available with the probe types and / or the outlet types shown below. If either the probe type or outlet type you require is not shown, please contact us to discuss your requirements.

PROBE TYPE



OUTLET TYPES



GAS SPECIFIC ADAPTORS



Twin Adaptor – BS Schrader – BS Probe

Gas	Part No.	Price
O ₂	673-9045-00	£120.73
MA4	673-9045-01	£120.73
MA7	673-9045-02	£120.73
O ₂ /N ₂ O	673-9045-03	£120.73
N ₂ O	673-9045-04	£120.73
Vacuum	673-9045-05	£120.73
CO ₂	673-9045-06	£120.73

- Twin Adaptor, 2x BS Schrader Outlets + BS Probe
- Designed for use with a 4 bar Oxygen supply
- Also available with different international outlet + probe types.



Twin Adaptor – BS Schrader – 3/8" Nut + Nipple

Gas	Part No.	Price
O ₂	673-9055-00	£109.73
MA4	673-9055-01	£109.73
MA7	673-9055-02	£109.73
O ₂ /N ₂ O	673-9055-03	£109.73
N ₂ O	673-9055-04	£109.73
Vacuum	673-9055-05	£109.73
CO ₂	673-9055-05	£109.73

- Twin Adaptor, 2x BS Schrader Outlets + 3/8" Nut + Nipple
- Designed for use with a 4 bar Oxygen supply
- Also available with different international outlet + probe types.

Each Gas Specific Adapter is available with the Probe Types and / or the Outlet Types shown below. If either the Probe Type or Outlet type you require is not shown, please contact us to discuss your requirements.

PROBE TYPE



OUTLET TYPES



HOSE ASSEMBLIES

LOW PRESSURE

**BS PROBE
BS SCHRADER
HOSE ASSEMBLY (3M)**



GAS	Part No.	Price
CO ₂	150-0002-00	£79.20
O ₂	151-0002-00	£79.20
O ₂ /N ₂ O	152-0002-00	£79.20
N ₂ O	153-0002-00	£79.20
MA7	154-0002-00	£79.20
MA4	155-0002-00	£79.20
VACUUM	156-0002-00	£79.20

**BS PROBE (90°)
BS SCHRADER
HOSE ASSEMBLY (3M)**



GAS	Part No.	Price
CO ₂	150-0102-00	£85.80
O ₂	151-0102-00	£85.80
O ₂ /N ₂ O	152-0102-00	£85.80
N ₂ O	153-0102-00	£85.80
MA7	154-0102-00	£85.80
MA4	155-0102-00	£85.80
VACUUM	156-0102-00	£85.80

**DISS NUT +TAIL
BS SCHRADER
HOSE ASSEMBLY (3M)**



GAS	Part No.	Price
CO ₂	150-1402-00	£79.20
O ₂	151-1402-00	£79.20
O ₂ /N ₂ O	152-1402-00	£79.20
N ₂ O	153-1402-00	£79.20
MA7	154-1402-00	£79.20
MA4	155-1402-00	£79.20
VACUUM	156-1402-00	£79.20

**DISS NUT + TAIL (90°)
BS SCHRADER
HOSE ASSEMBLY (3M)**



GAS	Part No.	Price
CO ₂	150-1502-00	£85.80
O ₂	151-1502-00	£85.80
O ₂ /N ₂ O	152-1502-00	£85.80
N ₂ O	153-1502-00	£85.80
MA7	154-1502-00	£85.80
MA4	155-1502-00	£85.80
VACUUM	156-1502-00	£85.80

**NIST PROBE + NUT
BS SCHRADER
HOSE ASSEMBLY (3M)**



GAS	Part No.	Price
CO ₂	150-0202-00	£79.20
O ₂	151-0202-00	£79.20
O ₂ /N ₂ O	152-0202-00	£79.20
N ₂ O	153-0202-00	£79.20
MA7	154-0202-00	£79.20
MA4	155-0202-00	£79.20
VACUUM	156-0202-00	£79.20

**NIST PROBE (90°)
BS SCHRADER
HOSE ASSEMBLY (3M)**



GAS	Part No.	Price
CO ₂	150-0302-00	£85.80
O ₂	151-0302-00	£85.80
O ₂ /N ₂ O	152-0302-00	£85.80
N ₂ O	153-0302-00	£85.80
MA7	154-0302-00	£85.80
MA4	155-0302-00	£85.80
VACUUM	156-0302-00	£85.80

HOSE ASSEMBLIES

LOW PRESSURE

BS PROBE DISS NUT + TAIL HOSE ASSEMBLY (3M)



GAS	Part No.	Price
CO ₂	150-0002-09	£47.30
O ₂	151-0002-09	£47.30
O ₂ /N ₂ O	152-0002-09	£47.30
N ₂ O	153-0002-09	£47.30
MA7	154-0002-09	£47.30
MA4	155-0002-09	£47.30
VACUUM	156-0002-09	£47.30

BS PROBE (90°) DISS NUT + TAIL HOSE ASSEMBLY (3M)



GAS	Part No.	Price
CO ₂	150-0102-09	£52.25
O ₂	151-0102-09	£52.25
O ₂ /N ₂ O	152-0102-09	£52.25
N ₂ O	153-0102-09	£52.25
MA7	154-0102-09	£52.25
MA4	155-0102-09	£52.25
VACUUM	156-0102-09	£52.25

BS PROBE DISS NUT + TAIL (90°) HOSE ASSEMBLY (3M)



GAS	Part No.	Price
CO ₂	150-0002-13	£52.25
O ₂	151-0002-13	£52.25
O ₂ /N ₂ O	152-0002-13	£52.25
N ₂ O	153-0002-13	£52.25
MA7	154-0002-13	£52.25
MA4	155-0002-13	£52.25
VACUUM	156-0002-13	£52.25

BS PROBE NIST PROBE + NUT HOSE ASSEMBLY (3M)



GAS	Part No.	Price
CO ₂	150-0002-01	£47.30
O ₂	151-0002-01	£47.30
O ₂ /N ₂ O	152-0002-01	£47.30
N ₂ O	153-0002-01	£47.30
MA7	154-0002-01	£47.30
MA4	155-0002-01	£47.30
VACUUM	156-0002-01	£47.30

BS PROBE (90°) NIST PROBE + NUT HOSE ASSEMBLY (3M)



GAS	Part No.	Price
CO ₂	150-0102-01	£52.25
O ₂	151-0102-01	£52.25
O ₂ /N ₂ O	152-0102-01	£52.25
N ₂ O	153-0102-01	£52.25
MA7	154-0102-01	£52.25
MA4	155-0102-01	£52.25
VACUUM	156-0102-01	£52.25

BS PROBE-NIST PROBE +NUT (90°) HOSE ASSEMBLY (3M)



GAS	Part No.	Price
CO ₂	150-0002-02	£52.25
O ₂	151-0002-02	£52.25
O ₂ /N ₂ O	152-0002-02	£52.25
N ₂ O	153-0002-02	£52.25
MA7	154-0002-02	£52.25
MA4	155-0002-02	£52.25
VACUUM	156-0002-02	£52.25

HOSE ASSEMBLIES

LOW PRESSURE

BS PROBE - 3/8" BSP NUT HOSE ASSEMBLY (3M)



GAS	Part No.	Price
CO2	150-0002-14	£46.75
O2	151-0002-14	£46.75
O2/N2O	152-0002-14	£46.75
N2O	153-0002-14	£46.75
MA7	154-0002-14	£46.75
MA4	155-0002-14	£46.75
VACUUM	156-0002-14	£46.75

BS PROBE (90°) 3/8" BSP NUT HOSE ASSEMBLY (3M)



GAS	Part No.	Price
CO2	150-0102-14	£52.25
O2	151-0102-14	£52.25
O2/N2O	152-0102-14	£52.25
N2O	153-0102-14	£52.25
MA7	154-0102-14	£52.25
MA4	155-0102-14	£52.25
VACUUM	156-0102-14	£52.25

BS PROBE 3/8" BSP NUT (90°) HOSE ASSEMBLY (3M)



GAS	Part No.	Price
CO2	150-0002-15	£52.25
O2	151-0002-15	£52.25
O2/N2O	152-0002-15	£52.25
N2O	153-0002-15	£52.25
MA7	154-0002-15	£52.25
MA4	155-0002-15	£52.25
VACUUM	156-0002-15	£52.25

BS PROBE – BS PROBE HOSE ASSEMBLY (3M)



GAS	Part No.	Price
CO ₂	150-0002-16	£46.75
O ₂	151-0002-16	£46.75
O ₂ /N ₂ O	152-0002-16	£46.75
N ₂ O	153-0002-16	£46.75
MA7	154-0002-16	£46.75
MA4	155-0002-16	£46.75
VACUUM	156-0002-16	£46.75

DISS NUT + TAIL- DISS NUT + TAIL HOSE ASSEMBLY (3M)



GAS	Part No.	Price
CO ₂	150-1402-09	£46.75
O ₂	151-1402-09	£46.75
O ₂ /N ₂ O	152-1402-09	£46.75
N ₂ O	153-1402-09	£46.75
MA7	154-1402-09	£46.75
MA4	155-1402-09	£46.75
VACUUM	156-1402-09	£46.75

NIST PROBE + NUT- NIST PROBE + NUT HOSE ASSEMBLY (3M)



GAS	Part No.	Price
CO ₂	150-0202-01	£46.75
O ₂	151-0202-01	£46.75
O ₂ /N ₂ O	152-0202-01	£46.75
N ₂ O	153-0202-01	£46.75
MA7	154-0202-01	£46.75
MA4	155-0202-01	£46.75
VACUUM	156-0202-01	£46.75

3/8" BSP NUT – 3/8" BSP NUT HOSE ASSEMBLY (3M)



GAS	Part No.	Price
CO ₂	150-1602-14	£46.75
O ₂	151-1602-14	£46.75
O ₂ /N ₂ O	152-1602-14	£46.75
N ₂ O	153-1602-14	£46.75
MA7	154-1602-14	£46.75
MA4	155-1602-14	£46.75
VACUUM	156-1602-14	£46.75

Create a Hose Assembly to suit your specific requirements by following the steps below to generate your **HOSE ASSEMBLY PART NUMBER**.


STEP 1 Choose the **GAS TYPE** that you require

STEP 2 Choose the **INLET / PROBE TYPE** you require

STEP 3 Choose the **HOSE LENGTH** you require

STEP 4 Choose the **OUTLET TYPE** you require

HOSE ASSEMBLY PART NUMBERING SYSTEM



GAS TYPES

O₂

O₂/N₂O

N₂O

MA7

MA4

VACUUM

CO₂

PROBE TYPES



BS Probe



BS Probe (90°)



NIST Probe + Nut



NIST Probe + Nut (90°)



AFNOR



AGA



Czech



SIS



DIN



UNIFOR



DISS HAND (F)



Ohmeda



Puritan



Chemetron



DISS Male



DISS NUT (F)



HOSE LENGTH

1M

2M

3M

4M

5M

6M

7M

OUTLET TYPES



BS Schrader



NIST Nut + Probe



NIST Probe + Nut (90°)



AFNOR



DIN



Czech



AGA



SIS



UNIFOR



DISS (M)




Ohmeda



Chemetron




Puritan



DISS HAND (F)



DISS NUT (F)





BS PROBE – TWIN BS SCHRADER ADAPTER – BAR HOSE ASSEMBLY

- 100mm bar Twin BS Schrader
- British Standard (other standards available upon request)
- 3M Hose

GAS	Part No.	Price
O ₂	673-9058-00	£157.30
O ₂ /N ₂ O	673-9059-00	£157.30
CO ₂	673-9088-00	£157.30
N ₂ O	673-9060-00	£157.30
MA7	673-9061-00	£157.30
MA4	673-9062-00	£157.30
VACUUM	673-9063-00	£157.30



BS PROBE – TWIN BS SCHRADER ADAPTOR – BAR HOSE ASSEMBLY + RAIL MOUNT

- 100mm bar Twin BS Schrader
- British Standard (other standards available upon request)
- 3M Hose
- Rail Mounted

GAS	Part No.	Price
O ₂	673-9064-00	£188.10
O ₂ /N ₂ O	673-9065-00	£188.10
CO ₂	673-9089-00	£188.10
N ₂ O	673-9066-00	£188.10
MA7	673-9067-00	£188.10
MA4	673-9068-00	£188.10
VACUUM	673-9069-00	£188.10



BS PROBE – TWIN BS SCHRADER ADAPTOR HOSE ASSEMBLY

- Twin BS Schrader Adaptor
- British Standard (other standards available upon request)
- 3M Hose

GAS	Part No.	Price
O ₂	673-9070-00	£163.90
O ₂ /N ₂ O	673-9071-00	£163.90
CO ₂	673-9090-00	£163.90
N ₂ O	673-9072-00	£163.90
MA7	673-9073-00	£163.90
MA4	673-9074-00	£163.90
VACUUM	673-9075-00	£163.90



BS PROBE – SINGLE BS SCHRADER + RAIL MOUNT HOSE ASSEMBLY

- Single BS Schrader
- British Standard (other standards available upon request)
- 3M Hose

GAS	Part No.	Price
O ₂	673-9076-00	£96.80
O ₂ /N ₂ O	673-9077-00	£96.80
CO ₂	673-9091-00	£96.80
N ₂ O	673-9078-00	£96.80
MA7	673-9079-00	£96.80
MA4	673-9080-00	£96.80
VACUUM	673-9081-00	£96.80



BS PROBE – TWIN BS SCHRADER + RAIL MOUNT HOSE ASSEMBLY

- Twin BS Schrader Adaptor
- British Standard (other standards available upon request)
- 3M Hose
- Rail Mounted

GAS	Part No.	Price
O ₂	673-9082-00	£165.00
O ₂ /N ₂ O	673-9083-00	£165.00
CO ₂	673-9092-00	£165.00
N ₂ O	673-9084-00	£165.00
MA7	673-9085-00	£165.00
MA4	673-9086-00	£165.00
VACUUM	673-9087-00	£165.00



Pin Index Yoke Filling Adaptor (1/4" BSP)

- Adaptor for filling Pin index cylinders
- Variety of thread types available – please specify
- Available for any medical gas type

GAS	Part No.	Price
O ₂	673-0621-00	£50.60
O ₂ /N ₂ O	673-0721-00	£50.60
N ₂ O	673-0821-00	£50.60
MA4	673-0921-00	£50.60
CO ₂	673-0021-00	£50.60



Pin Index Fixed Yoke Filling Adaptor (1/4" BSP)

- Adaptor for filling Pin index cylinders
- Variety of thread types available – please specify
- Available for any medical gas type

GAS	Part No.	Price
O ₂	673-1221-00	£50.60
O ₂ /N ₂ O	673-1321-00	£50.60
N ₂ O	673-1421-00	£50.60
MA4	673-1521-00	£50.60
CO ₂	673-1621-00	£50.60



Ferrule

- Stainless steel
- Different sizes available upon request

Part No.	Price
177-0723-00	£0.96



Heat Shrink Sleeving

- 25mm length
- Pack of 20

Colour	Part No.	Price
White	015-0014-00	£3.85
Black	015-0026-00	£3.85
Blue	015-0027-00	£3.85
Yellow	015-0028-00	£3.85
Clear	015-0029-00	£3.85
Grey	015-0030-00	£3.85



Medical Gas Hose

12.75mm OD x 6.35mm ID

- 30M coils
- Anti-static,
- Phthalate Free

GAS	Part No.	Price
O ₂	037-3101-03	£140.85
O ₂ /N ₂ O	037-3100-01	£140.85
CO ₂	037-3106-00	£140.85
N ₂ O	037-3300-00	£140.85
MA7	037-3500-00	£140.85
MA4	037-3101-02	£140.85
VACUUM	037-3013-00	£140.85



BS Probe

- ¼" Hose Barb (other sizes available)
- Stainless Steel
- Gas distinct 'collar'
- Meets British Standard
- MRI compatible

GAS	Part No.	Price
O ₂	177-0706-00	£15.40
MA4	177-0707-00	£15.40
MA7	177-1044-00	£15.40
N ₂ O	177-1043-00	£15.40
CO ₂	177-1140-00	£15.40
O ₂ /N ₂ O	177-0814-00	£15.40
Vacuum	177-1045-00	£15.40
Heliox	673-4006-00	£15.40



Right Angled BS Probe

- ¼" Hose Barb (other sizes available)
- Stainless Steel
- Gas distinct 'collar'
- Meets British Standard
- MRI compatible

GAS	Part No	Price
O ₂	673-4007-00	£29.70
MA4	673-4008-00	£29.70
MA7	673-4009-00	£29.70
N ₂ O	673-4010-00	£29.70
CO ₂	673-4220-00	£29.70
O ₂ /N ₂ O	673-4011-00	£29.70
Vacuum	673-4012-00	£29.70
Heliox	673-4013-00	£29.70



Direct BS Probe - 3/8" BSP Thread

- 3/8" BSP Thread
- Stainless Steel
- Gas distinct 'collar'
- Meets British Standard
- MRI compatible

GAS	Part No.	Price
O ₂	673-4014-00	£15.13
MA4	673-4015-00	£15.13
MA7	673-4016-00	£15.13
N ₂ O	673-4017-00	£15.13
CO ₂	673-4221-00	£15.13
O ₂ /N ₂ O	673-4018-00	£15.13
Vacuum	673-4019-00	£15.13
Heliox	673-4020-00	£15.13



Direct BS Probe M14x1 Thread

- M14x1 Thread
- Stainless Steel
- Gas distinct 'collar'
- Meets British Standard
- MRI compatible

GAS	Part No.	Price
O ₂	673-4021-00	£15.13
MA4	673-4022-00	£15.13
MA7	673-4023-00	£15.13
N ₂ O	673-4024-00	£15.13
CO ₂	673-4222-00	£15.13
O ₂ /N ₂ O	673-4025-00	£15.13



Direct BS Probe 9/16" UNF DISS

- 9/16" UNF DISS
- Stainless Steel
- Gas distinct 'collar'
- Meets British Standard
- MRI compatible

GAS	Part No	Price
O ₂	673-4026-00	£15.13
CO ₂	673-4223-00	£15.13

FITTINGS



NIST Nut + Probe Including O-Ring

- 1/4" Hose Barb (other sizes available)
- Stainless Steel (NIST PROBE)
- Gas specific engraving
- Gas specific indexing
- MRI compatible

GAS	Part No.	Price
O ₂	673-4027-00	£16.50
MA4	673-4028-00	£16.50
MA7	673-4029-00	£16.50
N ₂ O	673-4030-00	£16.50
CO ₂	673-4224-00	£16.50
O ₂ /N ₂ O	673-4031-00	£16.50
Vacuum	673-4032-00	£16.50
Heliox	673-4033-00	£16.50



NIST Probe Including O-Ring

- 1/4" Hose Barb (other sizes available)
- Stainless Steel
- Gas specific engraving
- Gas specific indexing
- MRI compatible

GAS	Part No.	Price
O ₂	177-0890-00	£9.35
MA4	177-1092-00	£9.35
MA7	673-4037-00	£9.35
N ₂ O	177-0893-00	£9.35
CO ₂	673-4225-00	£9.35
O ₂ /N ₂ O	673-4039-00	£9.35
Vacuum	673-4040-00	£9.35
Heliox	673-4041-00	£9.35



NIST Nut

- Plated brass nut
- Gas specific engraving on NIST Nut
- MRI compatible
- Stainless Steel

GAS	Part No.	Price
O ₂	177-0891-00	£8.25
MA4	177-1095-00	£8.25
MA7	673-4045-00	£8.25
N ₂ O	177-0896-00	£8.25
CO ₂	673-4226-00	£8.25
O ₂ /N ₂ O	673-4047-00	£8.25
Vacuum	673-4048-00	£8.25
Heliox	673-4049-00	£8.25



FEMALE NIST WITH HOSE TAIL

- Gas specific engraving
- MRI compatible
- Stainless Steel

GAS	Part No.	Price
O ₂	673-4051-00	£21.45
MA4	673-4052-00	£21.45
MA7	673-4053-00	£21.45
N ₂ O	673-4054-00	£21.45
CO ₂	673-4227-00	£21.45
O ₂ /N ₂ O	673-4055-00	£21.45
Vacuum	673-4056-00	£21.45



FEMALE NIST RIGHT ANGLED HOSE TAIL

- Gas specific Female NIST
- Gas specific engraving
- MRI compatible
- Stainless Steel

GAS	Part No.	Price
O ₂	673-4057-00	£33.55
MA4	673-4058-00	£33.55
MA7	673-4059-00	£33.55
N ₂ O	673-4060-00	£33.55
CO ₂	673-4228-00	£33.55
O ₂ /N ₂ O	673-4061-00	£33.55
Vacuum	673-4062-00	£33.55



FEMALE NIST 1/8" BSP THREAD

- Gas specific Female NIST
- Gas specific engraving
- MRI compatible
- Stainless Steel

GAS	Part No.	Price
O ₂	673-4063-00	£34.65
MA4	673-4064-00	£34.65
MA7	673-4065-00	£34.65
N ₂ O	673-4066-00	£34.65
CO ₂	673-4229-00	£34.64
O ₂ /N ₂ O	673-4067-00	£34.65
Vacuum	673-4068-00	£34.65

FITTINGS



AFNOR Probe

- French Standard
- Right Angle
- 1/4" Hose Barb (other sizes available)
- Stainless Steel

GAS	Part No.	Price
O ₂	673-4069-00	£33.00
MA4	673-4070-00	£33.00
N ₂ O	673-4071-00	£33.00
CO ₂	673-4230-00	£33.00
O ₂ /N ₂ O	673-4072-00	£33.00
Vacuum	673-4073-00	£33.00



AGA Probe

- Scandinavian Standard
- Right Angle
- 1/4" Hose Barb (other sizes available)
- Stainless Steel

GAS	Part No	Price
O ₂	673-4074-00	£33.00
MA4	673-4075-00	£33.00
N ₂ O	673-4076-00	£33.00
CO ₂	673-4231-00	£33.00
O ₂ /N ₂ O	673-4077-00	£33.00
Vacuum	673-4078-00	£33.00



Czech Probe

- Czech Republic Standard
- Right Angle
- 1/4" Hose Barb (other sizes available)
- Stainless Steel

GAS	Part No.	Price
O ₂	673-4079-00	£23.10
MA4	673-4080-00	£23.10
N ₂ O	673-4081-00	£23.10
CO ₂	673-4232-00	£23.10
O ₂ /N ₂ O	673-4082-00	£23.10
Vacuum	673-4083-00	£23.10



SIS Probe

- Australian Standard
- Right Angle
- 1/4" Hose Barb (other sizes available)
- Stainless Steel

GAS	Part No.	Price
O ₂	673-4084-00	£33.00
MA4	673-4085-00	£33.00
N ₂ O	673-4086-00	£33.00
CO ₂	673-4233-00	£33.00
O ₂ /N ₂ O	673-4087-00	£33.00
Vacuum	673-4088-00	£33.00



DIN Probe

- German Standard
- Right Angle
- 1/4" Hose Barb (other sizes available)
- Stainless Steel

GAS	Part No	Price
O ₂	673-4089-00	£23.10
MA4	673-4090-00	£23.10
N ₂ O	673-4091-00	£23.10
CO ₂	673-4234-00	£23.10
O ₂ /N ₂ O	673-4092-00	£23.10
Vacuum	673-4093-00	£23.10

FITTINGS



UNIFOR

- Italian Standard
- 1/4" Hose Barb (other sizes available)
- Stainless Steel

GAS	Part No.	Price
O ₂	673-4094-00	£33.00
MA4	673-4095-00	£33.00
N ₂ O	673-4096-00	£33.00
CO ₂	673-4235-00	£33.00
O ₂ /N ₂ O	673-4097-00	£33.00
Vacuum	673-4098-00	£33.00



DISS

- N. American Standard
- 1/4" Hose Barb (other sizes available)

GAS	Part No	Price
O ₂	673-4099-00	£16.50
MA4	673-4100-00	£16.50
N ₂ O	673-4101-00	£16.50
CO ₂	673-4236-00	£16.50
O ₂ /N ₂ O	673-4102-00	£16.50
Vacuum	673-4103-00	£16.50



Ohmeda

- N. American Standard
- 1/4" Hose Barb (other sizes available)

GAS	Part No.	Price
O ₂	673-4103-00	£16.50
MA4	673-4104-00	£16.50
N ₂ O	673-4105-00	£16.50
CO ₂	673-4237-00	£16.50
O ₂ /N ₂ O	673-4106-00	£16.50
Vacuum	673-4107-00	£16.50



Puritan

- N. American Standard
- 1/4" Hose Barb (other sizes available)

GAS	Part No.	Price
O ₂	673-4108-00	£33.00
MA4	673-4109-00	£33.00
N ₂ O	673-4110-00	£33.00
CO ₂	673-238-00	£33.00
O ₂ /N ₂ O	673-4111-00	£33.00
Vacuum	673-4112-00	£33.00



Chemetron

- N. American Standard
- 1/4" Hose Barb (other sizes available)

GAS	Part No	Price
O ₂	673-4113-00	£33.00
MA4	673-4114-00	£33.00
N ₂ O	673-4115-00	£33.00
CO ₂	673-239-00	£33.00
O ₂ /N ₂ O	673-4116-00	£33.00
Vacuum	673-4117-00	£33.00



DISS Male

- N. American Standard
- 1/4" Hose Barb (other sizes available)
- Stainless Steel

GAS	Part No.	Price
O ₂	673-4118-00	£12.37
MA4	673-4119-00	£12.37
N ₂ O	673-4120-00	£12.37
CO ₂	673-4240-00	£12.37
O ₂ /N ₂ O	673-4121-00	£12.37
Vacuum	673-4122-00	£12.37



**Single BS Schrader
1/4" BSPT Thread (Male)**

GAS	Part No.	Price
O ₂	673-4122-00	£30.25
MA4	673-4123-00	£30.25
MA7	673-4124-00	£30.25
N ₂ O	673-4125-00	£30.25
CO ₂	673-4241-00	£30.25
O ₂ /N ₂ O	673-4126-00	£30.25
Vacuum	673-4127-00	£30.25



**Single BS Schrader
1/8" BSPT Thread (Male)**

GAS	Part No.	Price
O ₂	673-4128-00	£30.25
MA4	673-4129-00	£30.25
MA7	673-4130-00	£30.25
N ₂ O	673-4131-00	£30.25
CO ₂	673-4242-00	£30.25
O ₂ /N ₂ O	673-4132-00	£30.25
Vacuum	673-4133-00	£30.25



**Single BS Schrader 1/4"
BSPT Thread Nut (Female)**

GAS	Part No.	Price
O ₂	673-4134-00	£30.25
MA4	673-4135-00	£30.25
MA7	673-4136-00	£30.25
N ₂ O	673-4137-00	£30.25
CO ₂	673-4243-00	£30.25
O ₂ /N ₂ O	673-4138-00	£30.25
Vacuum	673-4139-00	£30.25



**Single BS Schrader
3/8" Nut + Nipple**

GAS	Part No.	Price
O ₂	673-4140-00	£57.20
MA4	673-4141-00	£57.20
MA7	673-4142-00	£57.20
N ₂ O	673-4143-00	£57.20
CO ₂	673-4244-00	£57.20
O ₂ /N ₂ O	673-4144-00	£57.20
Vacuum	673-4145-00	£57.20



**Single BS Schrader +
Hose Tail 1/4"**

GAS	Part No.	Price
O ₂	673-4146-00	£42.35
MA4	673-4147-00	£42.35
MA7	673-4148-00	£42.35
N ₂ O	673-4149-00	£42.35
CO ₂	673-4244-00	£42.35
O ₂ /N ₂ O	673-4150-00	£42.35
Vacuum	673-4151-00	£42.35



**Single BS Schrader +
Hose Tail 1/8"**

GAS	Part No.	Price
O ₂	673-4152-00	£42.35
MA4	673-4153-00	£42.35
MA7	673-4154-00	£42.35
N ₂ O	673-4155-00	£42.35
CO ₂	673-4245-00	£42.35
O ₂ /N ₂ O	673-4156-00	£42.35
Vacuum	673-4157-00	£42.35

Each adaptor is available with the fittings below. If the fitting you require is not shown, please contact us to discuss your requirements.



1/4" BSPT THREAD
(MALE)



1/4" BSPT THREAD
(FEMALE)



HOSE TAIL



3/8" BSP NUT + NIPPLE

FITTINGS



Single AFNOR 1/4" BSPT Thread (Male)

GAS	Part No.	Price
O ₂	673-4160-00	£30.25
MA4	673-4161-00	£30.25
MA7	673-4162-00	£30.25
N ₂ O	673-4163-00	£30.25
CO ₂	673-4246-00	£30.25
O ₂ /N ₂ O	673-4164-00	£30.25
Vacuum	673-4165-00	£30.25



Single AGA 1/4" BSPT Thread (Male)

GAS	Part No.	Price
O ₂	673-4166-00	£30.25
MA4	673-4167-00	£30.25
MA7	673-4168-00	£30.25
N ₂ O	673-4169-00	£30.25
CO ₂	673-4247-00	£30.25
O ₂ /N ₂ O	673-4170-00	£30.25
Vacuum	673-4171-00	£30.25



Single CZECH 1/4" BSPT Thread (male)

GAS	Part No.	Price
O ₂	673-4172-00	£30.25
MA4	673-4173-00	£30.25
MA7	673-4174-00	£30.25
N ₂ O	673-4175-00	£30.25
CO ₂	673-4248-00	£30.25
O ₂ /N ₂ O	673-4176-00	£30.25
Vacuum	673-4177-00	£30.25



Single SIS 1/4" BSPT Thread (Male)

GAS	Part No.	Price
O ₂	673-4178-00	£30.25
MA4	673-4179-00	£30.25
MA7	673-4180-00	£30.25
N ₂ O	673-4181-00	£30.25
CO ₂	673-4249-00	£30.25
O ₂ /N ₂ O	673-4182-00	£30.25
Vacuum	673-4183-00	£30.25



Single DIN 1/4" BSPT Thread (Male)

GAS	Part No.	Price
O ₂	673-4184-00	£30.25
MA4	673-4185-00	£30.25
MA7	673-4186-00	£30.25
N ₂ O	673-4187-00	£30.25
CO ₂	673-4250-00	£30.25
O ₂ /N ₂ O	673-4188-00	£30.25
Vacuum	673-4189-00	£30.25



Single Unifor 1/4" BSPT Thread (Male)

GAS	Part No.	Price
O ₂	673-4190-00	£42.35
MA4	673-4191-00	£42.35
MA7	673-4192-00	£42.35
N ₂ O	673-4193-00	£42.35
CO ₂	673-4251-00	£42.35
O ₂ /N ₂ O	673-4194-00	£42.35
Vacuum	673-4195-00	£42.35

FITTINGS



Single DISS 1/4" BSPT Thread (Male)

GAS	Part No.	Price
O ₂	673-4196-00	£42.35
MA4	673-4197-00	£42.35
MA7	673-4198-00	£42.35
N ₂ O	673-4199-00	£42.35
CO ₂	673-4252-00	£42.35
O ₂ /N ₂ O	673-4200-00	£42.35
Vacuum	673-4201-00	£42.35



Single Ohmeda 1/4" BSPT Thread (Male)

GAS	Part No.	Price
O ₂	673-4202-00	£42.35
MA4	673-4203-00	£42.35
MA7	673-4204-00	£42.35
N ₂ O	673-4205-00	£42.35
CO ₂	673-4253-00	£42.35
O ₂ /N ₂ O	673-4206-00	£42.35
Vacuum	673-4207-00	£42.35



Single Chemetron 1/4" BSPT Thread (Male)

GAS	Part No.	Price
O ₂	673-4208-00	£42.35
MA4	673-4209-00	£42.35
MA7	673-4210-00	£42.35
N ₂ O	673-4211-00	£42.35
CO ₂	673-4254-00	£42.35
O ₂ /N ₂ O	673-4212-00	£42.35
Vacuum	673-4213-00	£42.35



Single Puritan 1/4" BSPT Thread (Male)

GAS	Part No.	Price
O ₂	673-4214-00	£42.35
MA4	673-4215-00	£42.35
MA7	673-4216-00	£42.35
N ₂ O	673-4217-00	£42.35
CO ₂	673-4255-00	£42.35
O ₂ /N ₂ O	673-4218-00	£42.35
Vacuum	673-4219-00	£42.35



BS Schrader Hose Tail

Size	Part No.	Price
1/4"	673-4158-00	£9.90
1/8"	673-4159-00	£9.90

COMPLETE VENTILATOR SYSTEM




Can be supplied with Huntleigh Vital Signs monitor which can be attached to the stand system.

- Complete ventilation System for use in Emergency and short / medium term stabilisation periods for patients in extremis.
- Gives a high level of ventilation control in a very compact and affordable package.
- Delivers 100% Oxygen or 50% concentrated which uses 2/3 less Oxygen.
- Ideal for field Hospitals and other Emergency applications.


Complete Package Part Number

670-0000-00 *Excludes Huntleigh Monitor

PRICE: POA



BEDSIDE MEDICAL OXYGEN SYSTEM



Oxygen Therapy System (Excluding Cylinder)

Part Number	Price
TBC	POA



MICROVENT™ WORLD PARAMEDIC RESUSCITATION KIT, WITH PIN INDEX AND BS SCHRADER FITTING

Part No. 670-0130-01

Price: £2,058.82

The Microvent™ World Resuscitation Kit is designed for professional emergency services use for when breathing has stopped, when breathing assistance is needed or if oxygen therapy is required.

Designed for easy use, suitable for adults + children.

Kit includes:

- 1x Microvent™ World Advanced Automatic + Manual Resuscitator With 2M Hose
- 1x Pin Index Regulator, 1x BS Schrader, Therapy (0-15LPM)
- 1x Size 5 Resuscitation Mask, 1x Size 4 Resuscitation Mask, 1x Size 3 Resuscitation Mask
- 1x Size 1, Size 2, Size 3 + Size 4 Oropharyngeal Airways
- 1x Catheter
- 1x Adult High Concentration Oxygen Therapy Mask
- 1x Child High Concentration Oxygen Therapy Mask



MICROVENT™ CPR PARAMEDIC RESUSCITATION KIT, WITH PIN INDEX AND BS SCHRADER FITTING

Part No. 670-0484-00

Price: £2,058.82

The Microvent™ CPR Resuscitation Kit is designed for professional emergency services use for when breathing has stopped, when breathing assistance is needed or if oxygen therapy is required.

Designed for easy use, suitable for adults + children.

Kit includes:

- 1x Microvent™ CPR Advanced Automatic + Manual Resuscitator With 2M Hose
- Set to 10 BPM with 1 Second Inspiratory Time
- 1x Pin Index Regulator, 1x BS Schrader, Therapy (0-15LPM)
- 1x Size 5 Resuscitation Mask, 1x Size 4 Resuscitation Mask, 1x Size 3 Resuscitation Mask
- 1x Size 1, Size 2, Size 3 + Size 4 Oropharyngeal Airways
- 1x Catheter
- 1x Adult High Concentration Oxygen Therapy Mask
- 1x Child High Concentration Oxygen Therapy Mask



OXYGEN THERAPY KIT AND MOUTH TO MASK RESUSCITATOR

Part No.	Price
670-0378-00	£344.85

Kit for giving Oxygen therapy or Oxygen enriched mouth to mask resuscitation. This kit includes:

- 1 x Pin Index therapy only Oxygen regulator
- 1 x High concentration Oxygen therapy masks (adult)
- 1 x High concentration Oxygen therapy masks (child)
- 4 x Oropharyngeal airways
- 1 x Pocket mask
- Supplied in soft, wipeable bag.

Oxygen Cylinder not included.



BAG VALVE MASK RESUSCITATION KIT I SUPPLIED WITH REGULATOR

Part No.	Price
670-0425-00	£453.75

- 1 x single use BVM with adult (5), small adult(4) and child(3) masks in silicone
- 1x Oxygen reservoir system
- 1 X Pocket Mask
- 2 x High concentration Oxygen therapy mask (adult)
- 2 x High concentration Oxygen therapy masks (child)
- 4 x Oropharyngeal airways
- Supplied in a soft, wipeable bag.

Oxygen Cylinder not included.



BAG VALVE MASK RESUSCITATION KIT II SUPPLIED WITHOUT REGULATOR

Part No.	Price
670-0425-90	£308.55

- 1 x single use BVM with adult (5), small adult (4), and child (3) masks in silicone
- 1 x Oxygen reservoir system
- 1 X Pocket Mask
- 2 x High concentration Oxygen therapy masks (adult)
- 2 x High concentration Oxygen therapy masks (child)
- 4 x Oropharyngeal airways
- Supplied in soft, wipeable bag.

Oxygen Cylinder not included.



MICROVENT™ INDUSTRIAL RESUSCITATION KIT (SOFT CASE).

Part No. 670-0309-00 **Price: £1,695.82**

MICROVENT™ INDUSTRIAL RESUSCITATION KIT (HARD CASE).

Part No. 670-0345-01 **Price: £1,724.85**

- Microvent™ advanced adult resuscitator with automatic and manual modes of operation
- 1 x Oxygen regulator, Pin Index, 1 x BS Schrader, therapy outlet
- 1 x Resuscitation mask (adult)
- 1 x High concentration Oxygen therapy mask (adult)
- 2 x Oropharyngeal airways
- 1 x Extension catheter mount (for connecting endo-tracheal tubes or other devices to the Microvent™ patient valve)
- Supplied in a soft, wipeable bag.

Oxygen Cylinder not included.



MICROVENT™ UTILITYVENT RESUSCITATION KIT.

Part No. 670-0340-01 **Price: £1,814.73**

The 'switch on and use' automatic resuscitator:

- Designed for ease of use where visibility and access may be reduced
- Pre-set automatic setting of 12 breaths per minute and 600ml breath size
- Patient assist valve for breathing synchronisation

The UtilityVENT™ Kit includes:

- 1 x UtilityVENT™ Resuscitator
- 1 x Oxygen regulator, Pin Index, 1 x BS Schrader
- 1 x Resuscitation mask (adult)
- 2 x Oropharyngeal airways
- Supplied in a soft, wipeable bag.

Oxygen Cylinder not included.



CONFINED SPACE RESCUE KIT - HARD CASE.

Part No. 670-0472-00
Price: £2,597.65

CONFINED SPACE RESCUE KIT -FABRIC CASE.

Part No. 670-0477-00
Price: £2,571.25

- 1 x Microvent™ advanced adult/child resuscitator with the 2m Oxygen supply hose
- 1 x Extension hose 10m (up to 30m available on request)
- Case for stowage and protection in adverse conditions
- 2 x High concentration Oxygen therapy mask (adult and child)
- 3 x Resuscitation masks (adult, small adult and child)
- 1 x Hook ring and head harness for attaching resuscitation mask
- 4 x Oropharyngeal airways
- 1 x Flexible ET tube mount
- 1 x Oxygen regulator, Pin Index, 1 x BS Schrader, therapy outlet
- 1 x Oxygen therapy Flow Selector (0-15 lpm)
- 1 x BA air regulator adapted to enable emergency use of breathing apparatus air with the Resuscitator.

Oxygen Cylinder not included.

Part No.
670-0498-00

Price
£1,695.82

MICROVENT™ SAILOR RESUSCITATION KIT, WITH PIN INDEX AND BS SCHRADER FITTING

The Microvent™ Sailor kit is designed for use in an emergency when breathing has stopped, when breathing assistance is needed or if Oxygen therapy is required.

Resuscitation Kit with:

- 1 x Microvent™ automatic and manual resuscitator with 2m Oxygen hose
- Designed for easy use
- Pre-set for adult casualties for ease of use in industry
- Manual setting ideal for 30:2 ventilation in CPR
- Automatic setting will provide consistent tireless ventilation
- Resuscitator is suitable for use in a non-breathable atmosphere
- Oxygen therapy through a therapy face mask or nasal cannula
- 1 x Oxygen regulator, Pin Index, 1 x BS Schrader, therapy outlet
- 2 x Oropharyngeal airways
- 1 x Catheter mount
- 1 x High concentration Oxygen therapy mask (adult)
- 1 x Resuscitation mask (adult)
- Supplied in soft, wipeable carry bag.

Oxygen Cylinder not included.



Part No.
670-0682-00

Price
£1,794.10

MICROVENT™ SAILOR II RESUSCITATION KIT, WITH PIN INDEX AND BS SCHRADER FITTING

The Microvent™ II Sailor kit is designed for use in an emergency when breathing has stopped, when breathing assistance is needed or if Oxygen therapy is required.

Resuscitation Kit with:

- 1 x Microvent™ automatic and manual resuscitator with 2m Oxygen hose
- Designed for easy use
- Pre-set for adult casualties for ease of use in industry
- Manual setting ideal for 30:2 ventilation in CPR
- Automatic setting will provide consistent tireless ventilation
- Oxygen therapy through a therapy face mask or nasal cannula
- 1 x Oxygen regulator, Pin Index, 1 x BS Schrader, therapy outlet
- 1 x Manual suction unit
- 1 x Pocket mask
- 1 x Cylinder holder for attaching Meditech 2L cylinder to other equipment such as a stretcher
- 4 x Oropharyngeal airways
- 1 x Catheter mount
- 1 x High concentration Oxygen therapy mask (adult)
- 1 x Resuscitation mask (adult)
- Supplied in soft, wipeable carry bag.

Oxygen Cylinder not included.





Part No. 900-3021-00 **Price: £551.10**

ENTONOX® DEMAND VALVE KIT WITH REGULATOR (2m HOSE)

- Standard Entonox® demand valve kit:
- Full length handle
 - 2m hose with BS Schrader probe
 - Entonox® demand valve kit bag
 - 5 x Mouth pieces/bacterial filters
 - 1 x Face mask (adult)
 - 1 x Face mask (child)
 - 1 x Entonox® regulator,
 - Customised kits are available. pin index, 1 x BS Schrader outlet

Part No. 900-1021-00 **Price: £387.75**

ENTONOX® DEMAND VALVE KIT WITHOUT REGULATOR (2m HOSE)

- Standard Entonox® demand valve kit:
- Full length handle
 - 2m hose with BS Schrader probe
 - No regulator
 - Entonox® demand valve kit bag
 - 1 x Face mask (adult)
 - 1 x Face mask (child)
 - Customised kits are available
 - 5 x Mouth pieces/bacterial filters.

Entonox Cylinder not included.



Part No. 910-3021-00 **Price: £551.10**

OXYGEN DEMAND VALVE KIT WITH REGULATOR (2m HOSE)

- Standard Oxygen Demand Valve Kit:
- Full length handle
 - 2m hose with BS Schrader probe
 - Oxygen demand valve kit bag
 - 5 x Mouth pieces/bacterial filters
 - 1 x Face mask (adult)
 - 1 x Face mask (child)
 - 1 x Oxygen regulator,
 - Customised kits are available. pin index, 1 x BS Schrader outlet

Part No. 910-0021-30 **Price: £366.03**

OXYGEN DEMAND VALVE KIT WITHOUT REGULATOR (2m HOSE)

- Standard Oxygen Demand Valve Kit:
- Full length handle
 - 2m hose with BS Schrader probe
 - No regulator
 - Oxygen demand valve kit bag
 - 1 x Face mask (adult)
 - 1 x Face mask (child)
 - Customised kits are available
 - 5 x Mouth pieces/bacterial filters.

Oxygen Cylinder not included.



MICROVENT™ MARINER KIT

Part No. 670-0413-00 Price: £2,837.45
MICROVENT™ MARINER RESCUE RESUSCITATOR WITH PIN INDEX AND BS SCHRADER FITTING

- 1 x Microvent™ advanced adult resuscitator with automatic and manual modes of operation. Fitted with 1m supply hose with BS Schrader probe.
- Automatic setting pre-set to Adult: 600ml breath at 12 breaths per minute
- 2 x Meditech Mariner Oxygen regulators (head end and foot end) with pin index inlet connector and BS Schrader outlet
- 1 x Pneumatic harness with switching shuttle valve fitted with BS medical Oxygen fittings
- 1 x Casualty harness with secured head harness, two cylinder storage pods with draw cords and zipped closure, resuscitator storage, waist strap
- 1 x Resuscitation mask (adult)
- 1 x Hook ring for attaching face mask to head harness
- Microvent™ Mariner carrying bag.

Oxygen Cylinder not included.



RESUSCITATION KIT FOR DENTAL PRACTITIONERS

A COMPREHENSIVE DENTAL RESUSCITATION KIT WHICH HAS BEEN DESIGNED TO MEET THE AIRWAY AND BREATHING REQUIREMENTS AS SET BY THE RESUSCITATION COUNCIL (UK) GUIDELINES.

Part No. 670-0684-00 Price: £659.45

Specific resuscitation equipment should be available immediately in all primary care dental premises. This equipment list is standardised throughout the UK. The resuscitation equipment in this dental kit is latex-free and for single-patient use. Replacement and spare parts can be found in the accessories pages.

- 1 x SL Pin Index therapy only regulator
- 1 x Pocket mask with Oxygen port
- 1 x Hand operated suction unit
- 5 x Oropharyngeal airways (size 0, 1, 2, 3, 4)
- 1 x BVM resuscitation system with size 4 mask (adult)
- 1 x Child BVM resuscitation system with size 3 mask
- 3 x Clear face mask for self-inflating bag (size 0, 1, 2)
- 1 x High concentration Oxygen therapy mask (adult)
- 1 x High concentration Oxygen therapy mask (child)
- 1 x Carry bag.

If you would like to customise your own kit, to suit your requirements, please contact our Sales Team to discuss further options.

Oxygen Cylinder not Included

BNOS Meditech Ltd. stock and supply a full range of Medical Oxygen Cylinders, ranging from 2L (D size) up to 40L (H size). They are available filled or empty with a Pin Index or Bullnose cylinder valve (other valves are available upon request). All Cylinders are filled to 200 bar and each Cylinder has a life span of approximately 3 years from the filling date. We have a Whole Sale Dealers license to distribute Medical Oxygen & Entonox.

Authorisation Holder Number: WDA (H) 44056



O₂ CYLINDER, 2L ALLOY, PIN INDEX VALVE

FILLED – Part No. 676-0362-00 Price: £242.00
EMPTY – Part No. 676-0361-00 Price: £183.00

Alloy Oxygen Cylinder supplied with medical Oxygen pin index valve:

- Filled to 200 bar (holding 400 litres of Oxygen)
- ISO 7866:1999 Gas Cylinders
- Empty Cylinder weight: 2.8kg
- Filled Weight (including Valve): 4kg
- Cylinder dimensions: Diameter – 111mm. Length – 420mm

O₂ CYLINDER, 2L ALLOY, BULLNOSE VALVE

FILLED – Part No. 676-0312-00 Price: £226.00
EMPTY – Part No. 676-0311-00 Price: £189.00

Alloy Oxygen Cylinder supplied with a medical Oxygen Bullnose valve:

- Filled to 200 bar (holding 400 litres of Oxygen)
- ISO 7866:1999 Gas Cylinders
- Empty Cylinder weight: 2.8kg
- Filled Weight (including Valve): 4kg
- Cylinder dimensions: Diameter – 111mm. Length – 420mm



O₂ CYLINDER, 5L STEEL, PIN INDEX VALVE

FILLED – Part No. 676-0025-00 Price: £251.00
EMPTY – Part No. 676-0026-00 Price: £160.00

Steel Oxygen Cylinder supplied with medical Oxygen pin index valve:

- Filled to 200 bar (holding 1000 litres of Oxygen)
- ISO 7866:1999 Gas Cylinders
- Empty Cylinder weight: 8.3kg
- Filled Weight (including Valve): 9.5kg
- Cylinder dimensions: Diameter – 140mm. Length – 560mm



O₂ CYLINDER, 10L STEEL, PIN INDEX VALVE

FILLED – Part No. 676-0029-00 Price: £285.00

EMPTY – Part No. 676-0030-00 Price: £200.00

Steel Oxygen Cylinder supplied with medical Oxygen pin index valve:

- Filled to 200 bar (holding 2000 litres of Oxygen)
- ISO 7866:1999 Gas Cylinders
- Empty Cylinder weight: 14.7kg
- Filled Weight (including Valve): 16kg
- Cylinder dimensions: Diameter – 140mm. Length – 800mm.

O₂ CYLINDER, 10L STEEL, BULLNOSE VALVE

FILLED – Part No. 676-0031-00 Price: £285.00

EMPTY – Part No. 676-0032-00 Price: £200.00

Steel Oxygen Cylinder supplied with a medical Oxygen Bullnose valve:

- Filled to 200 bar (holding 2000 litres of Oxygen)
- ISO 7866:1999 Gas Cylinders
- Empty Cylinder weight: 14.7kg
- Filled Weight (including Valve): 16kg
- Cylinder dimensions: Diameter – 140mm. Length – 800mm.



O₂ CYLINDER, 40L STEEL, PIN INDEX VALVE

FILLED – Part No. 676-0001-00 Price: £389.00

EMPTY – Part No. 676-0002-00 Price: £325.00

Steel Oxygen Cylinder supplied with medical Oxygen pin index valve:

- Filled to 200 bar (holding 8000 litres of Oxygen)
- ISO 7866:1999 Gas Cylinders
- Empty Cylinder weight: 36.2kg
- Filled Weight (including Valve): 40kg
- Cylinder dimensions: Diameter – 232mm. Length – 1300mm.

O₂ CYLINDER, 40L STEEL, BULLNOSE VALVE

FILLED – Part No. 676-0003-00 Price: £389.00

EMPTY – Part No. 676-0004-00 Price: £325.00

Steel Oxygen Cylinder supplied with a medical Oxygen Bullnose valve:

- Filled to 200 bar (holding 8000 litres of Oxygen)
- ISO 7866:1999 Gas Cylinders
- Empty Cylinder weight: 36.2kg
- Filled Weight (including Valve): 40kg
- Cylinder dimensions: Diameter – 232mm. Length – 1300mm.

Other sized available:
20L | 47L | 50L



ENTONOX CYLINDER, 2L ALLOY, PIN INDEX VALVE

FILLED – Part No. 676-0012-00 Price: £495.00

EMPTY – Part No. 676-0011-00 Price: £350.00

Alloy Entonox Cylinder supplied with Entonox Pin Index valve:

- Filled to 137 bar (holding 274 litres of Entonox)
- ISO 7866:1999 Gas Cylinders
- Empty Cylinder weight: 4kg
- Filled Weight (including Valve): 2.8kg
- Cylinder dimensions: Diameter – 111mm. Length – 420mm



ENTONOX CYLINDER, 10L ALLOY, PIN INDEX VALVE

FILLED – Part No. 676-0022-00 Price: £764.50

EMPTY – Part No. 676-0021-00 Price: £595.00

Alloy Entonox Cylinder supplied with medical Entonox Pin Index valve:

- Filled to 137 bar (holding 1370 litres of Entonox)
- ISO 7866:1999 Gas Cylinders
- Empty Cylinder weight: 14.7kg
- Filled Weight (including Valve): 16kg
- Cylinder dimensions: Diameter – 140mm. Length – 800mm.



CYLINDER RACK

Part No.	Price
674-0107-00	£91.74

- Secures Oxygen Cylinder to the wall – universal size.



CYLINDER SPINDLE KEY (PIN INDEX)

Part No.	Price
177-0112-00	£9.63



CYLINDER SPINDLE KEY (BULLNOSE)

Part No.	Price
674-0111-00	£8.80



BULLNOSE CYLINDER VALVE

17E THREAD		25E THREAD	
Part No.	Price	Part No.	Price
131-0145-00	£15.95	131-0185-00	£19.25



PIN INDEX CYLINDER VALVE

17E THREAD		25E THREAD	
Part No.	Price	Part No.	Price
131-0144-00	£15.95	131-0184-00	£19.25

ACCESSORIES

SINGLE USE OROPHARYNGEAL AIRWAYS



Size	Part No.	Price
Size 0.	673-0570-00	£1.43
Size 1.	673-0571-00	£1.43
Size 2.	673-0572-00	£1.43
Size 3.	673-0573-00	£1.43
Size 4.	673-0574-00	£1.43



WALL MOUNTING BRACKET

Part No.	Price
674-0133-00	£14.30

SINGLE USE BAG VALVE MASK



Adult Part No.	Price
673-0564-00	£25.91

Adult BVM resuscitation system with 1L bag with a pressure relief valve (40cm H₂O).

Child Part No.	Price
673-0565-00	£25.91

Child BVM resuscitation system with 550ml bag with a pressure relief valve (40cm H₂O).

- Both with 3m Oxygen tubing
- Oxygen reservoir bag with rapid bag recoil.

POCKET MASK



Part No.	Price
673-0569-00	£16.83

Pocket mask with supplementary O₂ Connector, head strap, filter & protective carry case. One way protection valve allows gas to pass from the rescuer to patient without any physical contact between the two, minimising the risk of cross contamination

HIGH CONCENTRATED THERAPY MASK



Adult Part No.	Price	Child Part No.	Price
673-9006-00	£8.47	673-9007-00	£8.47

- Oxygen therapy mask in PVC with PE connector
- Including 2m Oxygen tubing
- Elasticated adjustable head strap
- PE reservoir bag
- Single patient use.

ACCESSORIES



NASAL CANNULA

ADULT

Part No. Price
674-0086-00 £3.96

CHILD

Part No. Price
674-0113-00 £3.96

- Nasal cannula with straight prongs + tube (1.8M)



HOSE BARB CONVERTER (PLASTIC)

-Pack of 10

Part No. Price
674-0079-00 £20.68 (pack of 10)

Plastic converter, screws onto $\frac{9}{16}$ " thread and converts to fir tree style connector to accept plastic therapy hose.



THERAPY OUTLET FIR TREE – METAL

$\frac{9}{16}$ " UNF

Part No. Price
674-0137-00 £6.05

$\frac{3}{8}$ " BSP

Part No. Price
674-0103-00 £6.05

- Made from stainless steel
- Includes sealing O-rings
- Thread Type – $\frac{9}{16}$ " UNF or $\frac{3}{8}$ " BSP



THERAPY OUTLET ADAPTER

$\frac{9}{16}$ " UNF TO $\frac{3}{8}$ " BSP

Part No. Price
177-1029-03 £6.05

$\frac{9}{16}$ " UNF TO $\frac{1}{4}$ " BSP

Part No. Price
177-1040-00 £6.05

$\frac{3}{8}$ " BSP TO $\frac{9}{16}$ " UNF

Part No. Price
177-1053-00 £6.05

- Screw on adapter for threaded therapy outlets
- Made from stainless steel
- Includes sealing O-rings



BODOK SEAL

Part No. 033-1021-00 Price: £3.30

ACCESSORIES



HUMIDIFIER (single use or auto-clavable)

Single unit - single use

Part No.	Price
674-0047-00	£6.60

- Single-patient use
- Designed for long term durability
- With 350ml capacity

Carton of 60 (single use)

Part No.	Price
674-0047-60	£177.65

- Audible Safety pop alarm & a quad nut connector.
- Clear markings of maximum water level indicator

Single unit - auto-clavable | Part No. 674-0110-00 | Price: £21.95

- Designed for long term durability
- 350ml Capacity
- Audible Safety pop alarm & a quad nut connector.
- Clear markings of maximum water level indicator



BACTERIAL FILTER / MOUTH PIECE

Part No.	Price
673-9021-00	£8.75

Single use bacterial mouth pieces used to avoid cross contamination. To be used with our Entonox[®] and Oxygen demand valves.



MANUAL SUCTION UNIT

Part No. 673-0448-20	Price: £77.00
----------------------	---------------

- Manually operated suction unit
- Small and light for good portability
- Easy to use and direct
- Supplied with choice of catheters.

SUCTION JAR AND CATHETER - REPLACEMENTS

Part No. 674-0473-00	Price: £17.88
----------------------	---------------

Replacement suction jar with three different sizes of suction catheter.



RAIL CLAMP + WALL BRACKET

Part No.	Price
674-0126-00	£32.73

- Universal rail clamp
- Complete with wall bracket



RAIL MOUNTING BRACKET

Part No. 674-0129-00	Price: £14.30
----------------------	---------------

- 44MM Universal rail mounting bracket

ACCESSORIES



POLE CLAMP

Part No. 674-0134-00 Price: £54.73


- For use with 20-40mm diameter poles
- Complete with mounting bracket



RUBBER BOOT

Part No. 633-0027-00 Price: £5.75

- Protective rubber boot for Meditech Regulator Gauges



MEDICAL RAIL - 1M LENGTH

Part No. 674-0075-00 Price: £82.50

- Universal rail
- 3 x fixing screws & plugs and 2 end caps
- 1m long with 2m edge trim
- 1m cover strip.



OXYGEN SIGN

Part No. 674-0106-00 Price: £16.50



FABRIC BAG FOR ENTONOX KITS - WIPE CLEAN (SUPPLIED EMPTY)

Part No. 673-1018-00 Price £112.20

- Twin Carry Handles
- Reinforced, padded end panels
- Padded internal pockets to hold demand valve + accessories
- Clip on shoulder strap
- Marked 'ENTONOX KIT'
- Rubberised material covering base
- Padded sleeve to hold Entonox cylinder securely
- Double zip

Approximate dimensions: Diameter - 280mm. Length - 580mm.

The company Kolbrat & Bunz develop and produce RedVac and Medvac equipment for transport in rough terrain, for personal rescue & EMS use. All RedVac products are developed in close cooperation with Austrian and German Emergency Services, Alpine mountain rescue and other rescue organisations to ensure their equipment is manufactured to the highest quality and standard.



VACUUM MATTRESS

- Wide range of sizes and models available
- Multichamber system
- Wipe clean

Price: POA



VACUUM SPLINTS

- Wide range of sizes and models available
- Multichamber system
- Wipe clean

Price: POA



JUNIOR PEDIATRIC RESTRAINT SYSTEM

- For safe transport of children in Ambulances
- Compatible with most Ambulance gurneys
- Suitable for new-borns and children up to 7 years
- Max capacity: 130CM body length & 35KG

Price: POA

The company Kolbrat & Bunz develop and produce RedVac and Medvac equipment for transport in rough terrain, for personal rescue & EMS use. All RedVac products are developed in close cooperation with Austrian and German Emergency Services, Alpine mountain rescue and other rescue organisations to ensure their equipment is manufactured to the highest quality and standard.



FOOT PUMP - NON MAGNETIC

- Aluminium (also available in plastic)
- Length: 54cm
- Width: 6cm
- Weight: 1.4kg

Price: POA



HAND PUMP - NON MAGNETIC

- Aluminium (also available in plastic)
- Length: 30cm
- Width: 4cm
- Weight: 0.4kg

Price: POA



CARRY BAG

AVAILABLE FOR:

- Redvac junior
- Splint set
- Vacuum mattresses

Price: POA

Kohlbrat & Bunz develop and produce Medvac equipment for stanilising during MRI, CT examinations and operations.



ADULT MRI VACUUM MATTRESS

- Full body
- For safe immobilisation
- Easy to clean & disinfect
- Increased patient scans
- Patient is supported & imobilised
- Good insulation properties

Price: POA



CHILD 0-5 YEARS MRI VACUUM MATTRESS

- For safe immobilisation
- Easy to clean & disinfect
- Increased patient scans
- Patient is supported & imobilised
- Good insulation properties

Price: POA



INFANT 0-5 MONTHS MRI VACUUM MATTRESS

- For safe immobilisation
- Easy to clean & disinfect
- Increased patient scans
- Patient is supported & imobilised
- Good insulation properties

Price: POA

Kohlbrat & Bunz develop and produce Medvac equipment for stanilising during MRI, CT examinations and operations.



HAND & FOOT MRI VACUUM MATTRESS

- For safe immobilisation
- Easy to clean & disinfect
- Increased patient scans
- Patient is supported & imobilised
- Good insulation properties

Price: POA



ANKLE MRI VACUUM MATTRESS

- For safe immobilisation
- Easy to clean & disinfect
- Increased patient scans
- Patient is supported & imobilised
- Good insulation properties

Price: POA



KNEE, LEG & FOOT MRI VACUUM MATTRESS

- For safe immobilisation
- Easy to clean & disinfect
- Increased patient scans
- Patient is supported & imobilised
- Good insulation properties

Price: POA



HEAD IMMOBILISER

- For safe immobilisation
- Easy to clean & disinfect
- Increased patient scans
- Patient is supported & imobilised
- Good insulation properties

Price: POA

SERVICES AVAILABLE

PRODUCT SERVICING & MAINTENANCE

By customer request, regular servicing of all units is available. The units will be returned to Meditech where an initial examination takes place. Meditech will then send a quotation (including carriage) to the customer asking how they wish to proceed. Once the customer has given the go ahead, Meditech will proceed to carry out the check / service / rebuild of the unit. The unit will then be returned to the customer upon completion. We also offer a service contract set up which will be tailor made to suit the individual needs of the customer. A fixed-price calibration service on Microvent™ units is also available on request.

For more information on loan equipment and test kits please contact a member of our staff.

SERVICE TRAINING

We offer training courses so that a customer's engineer may learn to carry out various levels of service on their Meditech products. The basic levels are;

Level 1: Product familiarisation. Customer checks as detailed in user manual. Replacement of external parts without the use of tools.

Level 2: Calibration check using simple test equipment. Replacement of external parts including trigger and hose replacements. Use of simple tools will be required.

Level 3: Calibration checks using professional test equipment.

For more advanced levels of service training, please contact our sales team and they will be able to provide you with more information.

PLASTIC MOULDING

We have in-house capabilities to provide a wide range of Plastic Moulding Services to companies based Worldwide. With a dedicated team with over 50 years of experience, we are able to produce high quality plastic moulding products for a wide range of industries including; Medical, Engineering, Industrial & Aviation.

Plastic moulding Services we can offer;

- Sub contract plastic injection moulding services
- Two-part Liquid Silicone moulding
- Repair & maintenance of existing tools*
- Creation of new tools*

*Available in plastic and liquid silicone

PLATING

BNOS Meditech has a state of the art plating centre offering the best possible plating services to the Medical, Electronic, Marine, Aerospace, Engineering & manufacturing industries worldwide. We offer all types of electro-plating and polishing.

TECHNICAL SUPPORT

For any technical advice, support or pricing of our services please contact a member of staff via email or telephone and they will be able to help you.

TERMS AND CONDITIONS

1. DEFINITIONS AND INCORPORATION

- 1.1 In these Conditions, "Meditech" means B.N.O.S. Meditech Limited (Company Number 02368331) whose registered office is at Unit 9, Fifth Avenue, Bluebridge Industrial Estate, Halstead, Essex CO9 2SZ; "Contract" means any contract between Meditech and the Customer for the sale and purchase of the Goods and/or Services, including these Conditions; "Customer" means the party purchasing Goods and/or Services from Meditech; "Goods" means the medical equipment, kits and accessories supplied to the Customer by Meditech pursuant to the Contract; and "Services" means the repair, maintenance and related services be supplied to the Customer by Meditech pursuant to the Contract.
- 1.2 Subject to any variation pursuant to this Condition 1.2, the Contract will be on these Conditions to the exclusion of all other terms and conditions (including any terms or conditions which the Customer purports to apply under any purchase order, confirmation of order, specification or other document). No variation to these Conditions or any representations about the Goods and/or the Services shall have any effect unless expressly agreed in writing and signed by a duly authorised representative of Meditech.

2. ORDERS, DELIVERY, RISK AND TITLE

- 2.1 Each order for Goods and/or Services by the Customer from Meditech shall be deemed to be an offer by the Customer to purchase Goods and/or Services subject to these Conditions. No order placed by the Customer shall be deemed to be accepted by Meditech until a written acceptance of order is issued by Meditech or (if earlier) Meditech commences performance of the Contract.
- 2.2 The price of the Goods excludes carriage charges which the Customer shall pay in addition. Unless otherwise agreed in writing by Meditech, delivery of the Goods shall take place upon arrival of the Goods at the location specified in the Customer's order or otherwise agreed to in writing.
- 2.4 Any dates specified by Meditech for delivery of Goods and/or performance of the Services are an estimate only. Goods are at the risk of the Customer from the time of delivery but title in Goods shall not pass from Meditech until the Customer has paid all amounts due under the Contract. Until title in Goods passes to the Customer the Customer shall (a) hold Goods and each of them on a fiduciary basis as bailee for Meditech; (b) store Goods separately from all other goods in its possession (marked in such a way that they are clearly identified as Meditech's property); and (c) (subject to written agreement to the contrary by Meditech) not dispose of them. Further, until such time as title in Goods passes to the Customer, the Customer shall upon request forthwith deliver up such Goods to Meditech. The Customer hereby grants Meditech an irrevocable licence at any time to enter any premises where Goods are or may be stored in order to inspect them, and, until such time as title in Goods passes to the Customer, to recover them.

3. PRICE AND PAYMENT

- 3.1 The price for Goods and/or Services shall be the price set out in Meditech's quotation issued to the Customer, save that, unless otherwise stated in writing, any such quotation is valid for a period of 30 days from its date. Thereafter, a new quotation must be requested by the Customer.
- 3.2 All prices shall be exclusive of VAT, which if lawfully chargeable the Customer will pay in addition. Payment is due within [30] days of the date of Meditech's invoice. The Customer shall make all payments due under the Contract with out any deduction whether by way of set-off, counterclaim, discount, abatement or otherwise.
- 3.3 If the Customer fails to pay any sum due pursuant to the Contract the Customer will be liable to pay interest to Meditech on such sum from the due date for payment at the annual rate of [4]% above the base lending rate from time to time of NatWest Bank plc, accruing on a daily basis until payment is made, whether before or after any judgment.

4. WARRANTY AND LIABILITY

- 4.1 Subject to Condition 4.2 Meditech warrants that, provided the Customer and end-users follow Meditech's instructions as to the storage, installation, commissioning, use or maintenance of Goods or (if there are none) good industry practice, Goods will for a period of 12 months from the date of delivery, be free of material defects. If any Goods do not conform with such warranty Meditech shall at its option repair or replace such Goods as soon as reasonably practicable or refund the price of such Goods. This Condition 4.1 states the exclusive remedy of the Customer, and the entire liability of Meditech, for breach of such warranty.

TERMS AND CONDITIONS

- 4.2 Meditech shall not be liable for any breach of the above warranty given in relation to Goods unless the Customer gives written notice of the defect to Meditech within 30 days of the time when the Customer discovers or ought to have discovered the defect; and unless Customer (if asked to do so by Meditech) returns such Goods to Meditech at Customer cost for inspection by Meditech.
- 4.3 Meditech warrants that it will perform Services with reasonable care and skill. If it is shown to be in breach of such warranty in relation to particular Services, it shall at its own cost as soon as reasonably practicable re perform the relevant Services or instead (if it so decides) waive or refund to the Customer such of the fees payable or paid to Meditech by the Customer under the Contract as are attributable to those particular Services. This Condition 4.3 states the exclusive remedy of the Customer, and the entire liability of Meditech, for breach of such warranty.
- 4.4 Save as expressly set out herein, all warranties, conditions and other terms implied by statute or common law or otherwise are, to the fullest extent permitted by law, excluded from the Contract save that nothing in these Conditions excludes or limits the liability for fraud or for death or personal injury caused by negligence.
- 4.5 Subject to Condition 4.4, Meditech's total liability in contract, tort (including negligence or breach of statutory duty), misrepresentation or otherwise, arising out of or in connection with the Contract shall be limited to the price payable for the Goods and/or the Services which are the subject of the Customer's claim. Meditech shall not be liable to the Customer, in contract in tort or otherwise, for any loss of profits, business, revenue, goodwill or anticipated savings or for any indirect or consequential or economic loss whatsoever howsoever arising out of or in connection with the Contract.
- 4.6 Meditech shall not be liable or responsible for any loss or damage caused by delay in the performance or non-performance of its obligations hereunder where the same is occasioned by acts or omissions of the Customer, strikes, lock-out or trade dispute, difficulties in obtaining labour or materials or any other cause whatsoever beyond Meditech's reasonable control including but not limited to any act of God, war, terrorism, transport, civil disturbance, government or parliamentary restrictions, prohibitions or enactments of any kind. If the cause continues for more than 60 days either party may terminate the Contract by giving the other party no less than 30 days' written notice

5. GENERAL

- 5.1 The Customer may not assign or transfer its rights or delegate its duties under the Contract either in whole or in part. Meditech may assign the Contract in whole or part and/or delegate the performance of its obligations to third parties without reference to the Customer.
- 5.2 The failure of a party to exercise or enforce any right under the Contract shall not be deemed to be a waiver of that right nor operate to bar the exercise or enforcement of it at any time or times thereafter.
- 5.3 The Contract is not intended to confer a benefit on any third party, whether pursuant to the Contracts (Rights of Third Parties) Act 1999 or otherwise.
- 5.4 If any term of the Contract is held to be void or unenforceable by any reason of law, it shall be void or unenforceable to that extent only and no further and all other terms shall remain valid and fully enforceable.
- 5.5 Meditech may terminate the Contract forthwith by notice in writing to the Customer if the Customer ceases or threatens to cease to carry on its business or adopts a resolution for it to be wound up or if a petition is presented for the appointment of an administrator or a receiver or an administrative receiver is appointed in respect of any part of the Customer's undertaking or assets or if the Customer is unable to pay its debts within the meaning of Section 123 of the Insolvency Act 1986 (or any re-enactment of further enactment thereof). Termination of the Contract, for whatever cause, shall be without prejudice to the rights of either party accrued prior thereto, including without limitation any right to payment of any sum and any right to sue in respect of any antecedent breach of the Contract.
- 5.6 The Contract shall be governed by and construed in all respects in accordance with the laws of England and the parties submit to the exclusive jurisdiction of the English courts in respect of any dispute or claim arising out of or relating to the Contract.

medITECH

MEDICAL GAS EQUIPMENT

B.N.O.S. Meditech Ltd.

**9 Fifth Avenue, Bluebridge Industrial Estate,
Halstead, Essex CO9 2SZ, England, UK.**

Tel: +44 (0)1787 479475 • Fax: +44 (0)1787 477747

E-mail: sales@meditech.uk.com • www.meditech.uk.com



TO PURCHASE OR FOR FURTHER DETAILS:

**Contact us by email or telephone for more details,
including specific part numbers, prices, etc.**

All items in this range are sold with a return to manufacturer warranty against manufacturing defects. Specifications may be subject to change without notice. B.N.O.S. Meditech Ltd. terms and conditions of sale apply, available on application or from our website.

Distributor: